




Scan the code above to
learn more about our
Laparoscopic portfolio.

Laparoscopic Portfolio



symmetry surgical™

Laparoscopic



New Products	LAP-2
Accessories & Miscellaneous	LAP-22
Cables	LAP-31
Dissectors	LAP-32
MultiPak	LAP-32
ScopePak	LAP-32
Symmetry Surgical Laparoscopic Instruments	LAP-33
Symmetry Surgical RapidClean® Instruments	LAP-75
Symmetry Ultra™ Kleppinger Type Bipolar Forceps	LAP-90
Electrodes	LAP-92
Sterilization Solutions	LAP-94

Questions? Contact us at customerservice@symmetrysurgical.com
or learn more at symmetrysurgical.com



The Surgical Instrument Desk Reference

How to Use Our Catalog

The Symmetry Surgical catalog is designed with you in mind, to be your trusted resource for surgical instrumentation. We took great care in creating an easy-to-use, robust catalog and we think you'll be pleased with the design as you become more familiar with our broad portfolio.

Layout

The catalog is organized by surgical specialty **and** product function.

And we know that it is difficult to pick up an instrument and know what specialty it's listed in. To make your job easier, we've sorted the instruments in our portfolio by surgical function, too, including Energy, Laparoscopic, Retractors and Sterilization Solutions.

U.S. Customer Service

Your Symmetry Surgical Customer Service team wants to make ordering surgical instruments the easiest part of your day. We accept orders via phone, fax, email, our website, Global Healthcare Exchange (GHX) and through the mail.

PHONE

1.800.251.3000

(Available Monday through Friday, 8:30 AM–6:30 PM EST)

FAX

1.800.342.3272

EMAIL

customerservice@symmetrysurgical.com

ONLINE

symmetrysurgical.com

MAILING ADDRESS

Symmetry Surgical
Attn: Customer Service
3034 Owen Drive
Antioch, TN 37013

The Surgical Instrument Desk Reference

International Customer Service

PHONE

Symmetry Surgical Switzerland GmbH
+41 (0) 52 557 59 00
(Available Monday through Friday, 7:00-16:00 GMT)

Symmetry Surgical, Inc., USA
1.615.964.5672
(Available Monday through Friday, 8:30 AM-6:30 PM EST)

FAX

Symmetry Surgical Switzerland GmbH
+41 (0) 52 557 59 99

EMAIL

customerservice@symmetrysurgical.com
internationalorders@symmetrysurgical.com

ONLINE

symmetrysurgical.com

MAILING ADDRESS

Symmetry Surgical Switzerland GmbH
Bühlstrasse 27
CH-8200 Schaffhausen
Switzerland

Symmetry Surgical, Inc.
3034 Owen Drive
Antioch, TN 37013

INTERNATIONAL PAYMENT TERMS/RETURN GOODS POLICY/INSTRUMENT REPAIR & RESTORE PROGRAM/LITERATURE REQUESTS

Please contact your local distributor. Locate a local distributor contact, International Customer Service at +41 (0) 52 557 59 00 or 1.615.964.5672 (US).

Cardiovascular/Thoracic

Energy

Ear, Nose & Throat

General

Laparoscopic

Neurology

Obesity

OB/GYN

Ophthalmology

Orthopedics & Arthroscopy

Pediatrics

Rectal

Retractors

Spine

Urology

Sterilization Solutions

Alphabetical & Numeric Indexes

LAPAROSCOPIC

New Products

RapidClean® Glide Lock Instruments



RapidClean® Glide Lock Devices (Complete Device)

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
99-21302	5 mm Traumatic Grasper, 2x3 Teeth, Insert, Handle & Lumen, Standard Length (99-11302 & 99-00000)	5 mm	22 mm	37 mm
99-21302L	5 mm Traumatic Grasper, 2x3 Teeth, Insert, Handle & Lumen, Opti-Length® (99-11302L & 99-00000L)	5 mm	22 mm	49 mm
99-21323	5 mm Dorsey Fenestrated Grasper Insert, Handle & Lumen, Standard Length (99-11323 & 99-00000)	5 mm	30 mm	38 mm
99-21323L	5 mm Dorsey Fenestrated Grasper Insert, Handle & Lumen, Opti-Length® (99-11323L & 99-00000L)	5 mm	30 mm	49 mm
99-21333	5 mm Atraumatic Grasper Insert, Handle & Lumen, Standard Length (99-11333 & 99-00000)	5 mm	22 mm	37 mm
99-21333L	5 mm Atraumatic Grasper Insert, Handle & Lumen, Opti-Length® (99-11333L & 99-00000L)	5 mm	22 mm	49 mm
99-21345	5 mm Single Action Traumatic Grasper, 2x3 Teeth, Insert, Handle & Lumen, Standard Length (99-11345 & 99-00000)	5 mm	22 mm	37 mm
99-21345L	5 mm Single Action Traumatic Grasper, 2x3 Teeth, Insert, Handle & Lumen, Opti-Length® (99-11345L & 99-00000L)	5 mm	22 mm	48 mm
99-21346	5 mm Atraumatic Grasper Single Action Insert, Handle & Lumen, Standard Length (99-11346 & 99-00000)	5 mm	22 mm	37 mm
99-21346L	5 mm Atraumatic Grasper Single Action Insert, Handle & Lumen, Opti-Length® (99-11346L & 99-00000L)	5 mm	22 mm	48 mm
99-21364	5 mm Tenaculum Insert, Handle & Lumen, Standard Length (99-11364 & 99-00000)	5 mm	23 mm	37 mm
99-21364L	5 mm Tenaculum Insert, Handle & Lumen, Opti-Length® (99-11364L & 99-00000L)	5 mm	23 mm	49 mm
99-21369	5 mm Fenestrated Bowel Grasper Insert, Handle & Lumen, Standard Length (99-11369 & 99-00000)	5 mm	20 mm	37 mm
99-21369L	5 mm Fenestrated Bowel Grasper Insert, Handle & Lumen, Opti-Length® (99-11369L & 99-00000L)	5 mm	20 mm	48 mm

New Products

RapidClean® Glide Lock Instruments

Replacement Handle & Lumen Assemblies for RapidClean® Glide Lock

CODE	DESCRIPTION	DIAMETER	TIP SIZE
99-00000	5 mm Take-Apart Glide Lock Handle & Tube Assembly, Standard Length	5 mm	28 mm
99-00000L	5 mm Take-Apart Glide Lock Handle & Tube Assembly, Opti-Length®	5 mm	28 mm



99-11302

Replacement Insert Assemblies for RapidClean® Glide Lock

CODE	DESCRIPTION	DIAMETER	TIP SIZE
99-11302	5 mm Traumatic Grasper, 2x3 Teeth, Insert, Standard Length	5 mm	22 mm
99-11302L	5 mm Traumatic Grasper, 2x3 Teeth, Insert, Opti-Length®	5 mm	22 mm
99-11323	5 mm Dorsey Fenestrated Grasper Insert, Standard Length	5 mm	30 mm
99-11323L	5 mm Dorsey Fenestrated Grasper Insert, Opti-Length®	5 mm	30 mm
99-11333	5 mm Atraumatic Grasper Insert, Standard Length	5 mm	22 mm
99-11333L	5 mm Atraumatic Grasper Insert, Opti-Length®	5 mm	22 mm
99-11345	5 mm Single Action Traumatic Grasper, 2x3 Teeth, Insert, Standard Length	5 mm	22 mm
99-11345L	5 mm Single Action Traumatic Grasper, 2x3 Teeth, Insert, Opti-Length®	5 mm	22 mm
99-11346	5 mm Atraumatic Grasper Single Action Insert, Standard Length	5 mm	22 mm
99-11346L	5 mm Atraumatic Grasper Single Action Insert, Opti-Length®	5 mm	22 mm
99-11364	5 mm Tenaculum Insert, Standard Length	5 mm	23 mm
99-11364L	5 mm Tenaculum Insert, Opti-Length®	5 mm	23 mm
99-11369	5 mm Fenestrated Bowel Grasper Insert, Standard Length	5 mm	20 mm
99-11369	5 mm Fenestrated Bowel Grasper Insert, Opti-Length®	5 mm	20 mm

New Products

Symmetry Surgical Laparoscopic Instruments

Single Piece Laparoscopic Instruments are available with multiple handle configurations.



5mm Allis Forcep

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1111	Monopolar Handle	5 mm	12 mm	33 mm
98-1112	Monopolar Handle with Rotation	5 mm	12 mm	35 mm
98-1113	Ratcheting Handle	5 mm	12 mm	33 mm
98-1114	Ratcheting Handle with Rotation	5 mm	12 mm	35 mm
98-114L	Ractheting Handle with Rotation Opti-Length®, 45	5 mm	12 mm	45 mm

New Products

Symmetry Surgical Laparoscopic Instruments



5mm Long Allis Forcep

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1115	Monopolar Handle with Rotation	5 mm	28 mm	37 mm
98-1116	Ratcheting Handle with Rotation	5 mm	28 mm	37 mm
98-1117	Slide Lock Handle with Rotation	5 mm	28 mm	37 mm
98-1118	Standard Handle with Rotation	5 mm	28 mm	37 mm



5mm Traumatic Grasper with 2x3 Teeth

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1119	Auto Closure Handle with Rotation	5 mm	22 mm	36 mm
98-1120	Ratcheting Handle with Rotation	5 mm	22 mm	36 mm
95-1048SL-DA	Slide Lock Handle	5 mm	22 mm	36 mm
95-1048SL-DA-ROT	Slide Lock Handle with Rotation	5 mm	22 mm	36 mm
98-1121L	Slide Lock Handle with Rotation Opti-Length®, 45	5 mm	22 mm	45 mm
98-1122	Standard Handle with Rotation	5 mm	22 mm	36 mm



5mm Glassman Forcep

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1123	Monopolar Handle with Rotation	5 mm	40 mm	38 mm
98-1124	Ratcheting Handle	5 mm	40 mm	36 mm
98-1125	Ratcheting Handle with Rotation	5 mm	40 mm	38 mm
98-1125L	Ratcheting Handle with Rotation Opti-Length®, 45	5 mm	40 mm	47 mm
98-1126	Standard Handle	5 mm	40 mm	36 mm
98-1127	Slide Lock Handle with Rotation	5 mm	40 mm	38 mm
98-1127L	Slide Lock Handle with Rotation Opti-Length®, 45	5 mm	40 mm	47 mm
98-1128	Standard Handle with Rotation	5 mm	40 mm	38 mm
98-1128L	Standard Handle with Rotation Opti-Length®, 45	5 mm	40 mm	47 mm

New Products

Symmetry Surgical Laparoscopic Instruments



5mm Mckernan Tri-Grasping Fenestrated Forcep

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1129	Ratcheting Handle	5 mm	21 mm	34 mm
98-1130	Ratcheting Handle with Rotation	5 mm	21 mm	36 mm
98-1130L	Ratcheting Handle with Rotation Opti-Length®, 45	5 mm	21 mm	45 mm
98-1131	Slide Lock Handle with Rotation	5 mm	21 mm	36 mm
98-1132	Standard Handle with Rotation	5 mm	21 mm	36 mm



5mm Alligator Forcep

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1133	Monopolar Handle	5 mm	13 mm	36 mm
90-1072	Ratcheting Handle	5 mm	13 mm	36 mm
90-1072L	Monopolar Handle Opti-Length®, 45	5 mm	13 mm	44 mm
98-1134	Ratcheting Handle with Rotation	5 mm	13 mm	35 mm



5mm Pointed Dissector with Cup

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
90-1033	Monopolar Handle	5 mm	12 mm	35 mm
90-1033E	Monopolar Handle Opti-Length®, 45	5 mm	12 mm	44 mm
98-1135	Monopolar Handle with Rotation	5 mm	12 mm	35 mm
98-1136	Standard Handle with Rotation	5 mm	12 mm	35 mm



5mm Curved Dissecting Forcep

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1137	Monopolar Handle	5 mm	19 mm	36 mm
98-1138	Monopolar Handle with Rotation	5 mm	19 mm	36 mm
98-1139	Standard Handle with Rotation	5 mm	19 mm	36 mm

New Products

Symmetry Surgical Laparoscopic Instruments



5mm Short Curved Dissector

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1140	Monopolar Handle	5 mm	13 mm	33 mm
98-1141	Monopolar Handle with Rotation	5 mm	13 mm	35 mm



5mm Mixer Dissector

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1142	Monopolar Handle	5 mm	19 mm	33 mm
98-1143	Monopolar Handle with Rotation	5 mm	19 mm	36 mm
98-1144	Standard Handle with Rotation	5 mm	19 mm	36 mm



5mm Mini Right Angle Dissector

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1145	Auto Closure Handle with Rotation	5 mm	12 mm	35 mm
98-1146	Monopolar Handle with Rotation	5 mm	12 mm	35 mm
98-1146L	Monopolar Handle with Rotation Opti-Length®, 45	5 mm	12 mm	44 mm
98-1147	Ratcheting Handle with Rotation	5 mm	12 mm	35 mm
98-1148	Slide Lock Handle	5 mm	12 mm	35 mm
98-1149	Slide Lock Handle with Rotation	5 mm	12 mm	35 mm
98-1150	Standard Handle with Rotation	5 mm	12 mm	35 mm
98-1150L	Standard Handle with Rotation Opti-Length®, 45	5 mm	12 mm	44 mm



5mm Micro Dolphin Dissecting Forcep

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1151	Monopolar Handle	5 mm	11 mm	33 mm
98-1152	Monopolar Handle with Rotation	5 mm	11 mm	35 mm

New Products

Symmetry Surgical Laparoscopic Instruments



5mm Claw Forcep with 2x3 Teeth

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
90-1005	Auto Closure Handle	5 mm	22 mm	34 mm
98-1153	Ratcheting Handle	5 mm	22 mm	34 mm
98-1154	Ratcheting Handle with Rotation	5 mm	22 mm	36 mm
98-1154L	Slide Lock Handle with Rotation, 45	5 mm	22 mm	36 mm
98-1155	Slide Lock Handle with Rotation	5 mm	22 mm	36 mm



5mm Biopsy Forcep with Two Teeth

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1156	Auto Closure Handle	5 mm	11 mm	33 mm
90-1071	Monopolar Handle	5 mm	11 mm	33 mm
90-1071L	Monopolar Handle Opti-Length®, 45	5 mm	11 mm	44 mm



5mm Dorsey Fenestrated Grasper

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
97-1031RC	Ratcheting Handle with Rotation	5 mm	30 mm	37 mm
97-1031RCL	Ratcheting Handle with Rotation Opti-Length®, 45	5 mm	30 mm	46 mm
95-1031SL	Slide Lock Handle Opti-Length®, 45	5 mm	30 mm	46 mm
98-1157	Slide Lock Handle with Rotation	5 mm	30 mm	37 mm
98-1157L	Slide Lock with Rotation Opti-Length®, 45	5 mm	30 mm	46 mm
98-1158	Standard Handle with Rotation	5 mm	30 mm	37 mm

New Products

Symmetry Surgical Laparoscopic Instruments



5mm Endo Babcock Forcep

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1159	Monopolar Handle with Rotation	5 mm	29 mm	37 mm
97-1089RC	Ratcheting Handle with Rotation	5 mm	29 mm	37 mm
98-1160	Slide Lock Handle with Rotation	5 mm	29 mm	37 mm
98-1161	Standard Handle with Rotation	5 mm	29 mm	37 mm



5mm DeBakey Tissue Forcep

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1162	Monopolar Handle with Rotation	5 mm	38 mm	38 mm
98-1163L	Ratcheting Handle Opti-Length®, 45	5 mm	38 mm	45 mm
98-1164	Ratcheting Handle with Rotation	5 mm	38 mm	38 mm
98-1164L	Ratcheting Handle with Rotation Opti-Length®, 45	5 mm	38 mm	45 mm
98-1165	Standard Handle	5 mm	38 mm	35 mm
98-1165L	Standard Handle Opti-Length®, 45	5 mm	38 mm	45 mm
98-1166	Slide Lock Handle with Rotation	5 mm	38 mm	38 mm
98-1167L	Standard Handle with Rotation Opti-Length®, 45	5 mm	38 mm	45 mm



5mm Grasping Forcep

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1168	Auto Closure Handle with Rotation	5 mm	13 mm	35 mm
98-1169	Monopolar Handle	5 mm	13 mm	33 mm
98-1170	Monopolar Handle with Rotation	5 mm	13 mm	35 mm
98-1170L	Monopolar Handle with Rotation, 45	5 mm	13 mm	35 mm
98-1171	Ratcheting Handle	5 mm	13 mm	33 mm
98-1172	Ratcheting Handle with Rotation	5 mm	13 mm	35 mm
98-1172L	Ratcheting Handle with Rotation, 45	5 mm	13 mm	35 mm
98-1173	Standard Handle	5 mm	13 mm	33 mm
98-1174	Slide Lock Handle with Rotation	5 mm	13 mm	35 mm
98-1175	Standard Handle with Rotation	5 mm	13 mm	35 mm

New Products

Symmetry Surgical Laparoscopic Instruments



5mm Fenestrated Grasping Forcep

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1176	Monopolar Handle with Rotation	5 mm	22 mm	36 mm
98-1177	Ratcheting Handle with Rotation	5 mm	22 mm	36 mm
98-1178	Slide Lock Handle with Rotation	5 mm	22 mm	36 mm



5mm Micro Dissecting Forcep

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
90-1044	Monopolar Handle	5 mm	13 mm	33 mm
90-1044L	Monopolar Handle Opti-Length®, 45	5 mm	13 mm	45 mm
97-1042RC	Monopolar Handle with Rotation	5 mm	13 mm	35 mm
98-1179	Ratcheting Handle with Rotation	5 mm	13 mm	35 mm



5mm Babcock Grasping Forcep

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1180	Monopolar Handle with Rotation	5 mm	25 mm	37 mm
98-1181	Ratcheting Handle	5 mm	25 mm	34 mm
98-1182	Ratcheting Handle with Rotation	5 mm	25 mm	37 mm
98-1183	Slide Lock Handle with Rotation	5 mm	25 mm	37 mm



5mm Dolphin Nose Dissecting Forcep with Cup

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
90-1030	Monopolar Handle	5 mm	16 mm	36 mm
90-1030L	Monopolar Handle Opti-Length®, 45	5 mm	16 mm	45 mm
97-1030RC	Monopolar Handle with Rotation	5 mm	16 mm	36 mm


New Products

Symmetry Surgical Laparoscopic Instruments



5mm Cobra Jaw Grasping Forcep

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1184	Ratcheting Handle	5 mm	16 mm	33 mm
90-1069	Ratcheting Handle	5 mm	13 mm	35 mm
90-1069RL	Ratcheting Handle with Rotation Opti-Length®, 45	5 mm	13 mm	45 mm



5mm Atraumatic Grasper (Double Action)

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1185	Auto Closure Handle with Rotation	5 mm	22 mm	36 mm
98-1186	Monopolar Handle with Rotation	5 mm	22 mm	36 mm
98-1186L	Monopolar Handle with Rotation Opti-Length®, 45	5 mm	22 mm	45 mm
98-1187	Ratcheting Handle	5 mm	22 mm	34 mm
98-1188	Ratcheting Handle with Rotation	5 mm	22 mm	36 mm
98-1188L	Ratcheting Handle with Rotation Opti-Length®, 45	5 mm	22 mm	45 mm
98-1189	Slide Lock Handle	5 mm	22 mm	34 mm
95-1049SL-DA-45	Slide Lock Handle Opti-Length®, 45	5 mm	22 mm	46 mm
95-1049SL-DA-ROT	Slide Lock Handle with Rotation	5 mm	22 mm	36 mm
98-1109L	Slide Lock Handle with Rotation Opti-Length®, 45	5 mm	22 mm	45 mm
98-1191	Standard Handle with Rotation	5 mm	22 mm	36 mm
98-1191L	Standard Handle with Rotation Opti-Length®, 45	5 mm	22 mm	45 mm



5mm Paddle Babcock

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1192	Ratcheting Handle with Rotation	5 mm	29 mm	37 mm
98-1193	Standard Handle with Rotation	5 mm	29 mm	37 mm
98-1194	Monopolar Handle with Rotation	5 mm	29 mm	37 mm
98-1195	Slide Lock Handle with Rotation	5 mm	29 mm	37 mm

New Products

Symmetry Surgical Laparoscopic Instruments



5mm Tri-Grasper

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
90-1050	Ratcheting Handle	5 mm	16 mm	33 mm
90-1050L	Ratcheting Handle Opti-Length®, 45	5 mm	16 mm	45 mm
97-1050RC	Ratcheting Handle with Rotation	5 mm	16 mm	36 mm
95-1050SL-DA-ROT	Slide Lock Handle with Rotation	5 mm	16 mm	36 mm
95-1050SL-DA-ROT-45	Slide Lock with Rotation Opti-Length®, 45	5 mm	16 mm	45 mm
98-1196	Standard Handle with Rotation	5 mm	16 mm	36 mm



5mm Long Grasping Forcep

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1197	Auto Closure Handle with Rotation	5 mm	19 mm	36 mm
98-1198	Monopolar Handle with Rotation	5 mm	19 mm	36 mm
98-1198L	Monopolar Handle with Rotation, 45	5 mm	19 mm	36 mm
98-1199	Ratcheting Handle with Rotation	5 mm	19 mm	36 mm
98-1199L	Ratcheting Handle with Rotation, 45	5 mm	19 mm	36 mm
98-1200	Slide Lock Handle with Rotation	5 mm	19 mm	36 mm
98-1201	Standard Handle with Rotation	5 mm	19 mm	36 mm



5mm Maxi Fenestrated Grasper

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1202	Monopolar Handle with Rotation	5 mm	18 mm	36 mm
98-1203	Ratcheting Handle	5 mm	18 mm	35 mm
98-1204	Ratcheting Handle with Rotation	5 mm	18 mm	36 mm
98-1204L	Ratcheting Handle with Rotation Opti-Length®, 45	5 mm	18 mm	45 mm
98-1205	Slide Lock Handle with Rotation	5 mm	18 mm	36 mm

New Products

Symmetry Surgical Laparoscopic Instruments



5mm Retrograde Grasping Forcep

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1206	Ratcheting Handle	5 mm	13 mm	33 mm
98-1207	Ratcheting Handle with Rotation	5 mm	13 mm	35 mm
98-1208	Slide Lock Handle with Rotation	5 mm	13 mm	35 mm



5mm Kocher Grasping Forcep

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1209	Auto Closure Handle with Rotation	5 mm	13 mm	35 mm
98-1210	Ratcheting Handle with Rotation	5 mm	13 mm	35 mm
98-1211	Ratcheting Handle	5 mm	13 mm	35 mm
98-1212	Standard Handle with Rotation	5 mm	13 mm	35 mm



5mm Curved Maryland Dissector

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
90-1067	Monopolar Handle	5 mm	23 mm	34 mm
90-1067-L	Monopolar Handle Opti-Length®, 45	5 mm	23 mm	45 mm
97-1067RC	Monopolar Handle with Rotation	5 mm	23 mm	36 mm
98-1286L	Monopolar Handle with Rotation, 45	5 mm	23 mm	36 mm
95-1067SL-ROT-45	Slide Lock Handle with Rotation Opti-Length®, 45	5 mm	23 mm	45 mm
98-1200	Slide Lock Handle with Rotation	5 mm	23 mm	36 mm

New Products

Symmetry Surgical Laparoscopic Instruments



5mm Dolphin Nose Dissecting Forcep with Cup

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1213	Monopolar Handle with Rotation	5 mm	13 mm	35 mm
98-1214	Ratcheting Handle with Rotation	5 mm	13 mm	35 mm
98-1215	Slide Lock Handle with Rotation	5 mm	13 mm	35 mm
98-1216	Standard Handle with Rotation	5 mm	13 mm	35 mm



5mm Single Action Traumatic Grasper

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
97-1056RC	Ratcheting Handle with Rotation	5 mm	22 mm	36 mm
95-1048SL	Slide Lock Handle	5 mm	22 mm	36 mm
95-1048SL-45	Slide Lock Handle Opti-Length®, 45	5 mm	22 mm	45 mm
98-1217	Slide Lock Handle with Rotation	5 mm	22 mm	36 mm
98-1217L	Slide Lock Handle with Rotation Opti-Length®, 45	5 mm	22 mm	45 mm



5mm Single Action Atraumatic Grasper

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
97-1051RC	Ratcheting Handle with Rotation	5 mm	22 mm	36 mm
98-1218L	Ratcheting Handle with Rotation Opti-Length®, 45	5 mm	22 mm	45 mm
95-1049SL	Slide Lock Handle	5 mm	22 mm	37 mm
95-1049SL-45	Slide Lock Handle Opti-Length®, 45	5 mm	22 mm	45 mm
95-1049SL-ROT	Slide Lock Handle with Rotation	5 mm	22 mm	36 mm
95-1049SL-ROT-45	Slide Lock Handle with Rotation Opti-Length®, 45	5 mm	22 mm	45 mm
98-1220	Standard Handle with Rotation	5 mm	22mm	36 mm

New Products

Symmetry Surgical Laparoscopic Instruments



5mm Grasping Forcep with Cup

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1221	Auto Closure Handle with Rotation	5 mm	12 mm	35 mm
90-1035	Monopolar Handle	5 mm	12 mm	33 mm
90-1035L	Monopolar Handle OptiLength, 45	5 mm	12 mm	45 mm
97-1035RC	Monopolar Handle with Rotation	5 mm	12 mm	35 mm
98-1222L	Monopolar Handle with Rotation Opti-Length®, 45	5 mm	12 mm	44 mm
97-1055RC	Ratcheting Handle with Rotation	5 mm	12 mm	35 mm
98-1223L	Ratcheting Handle with Rotation Opti-Length®	5 mm	12 mm	44 mm
98-1224	Slide Lock Handle with Rotation	5 mm	12 mm	35 mm
98-1225	Standard Handle with Rotation	5 mm	12 mm	35 mm
98-1225L	Standard Handle with Rotation Opti-Length®, 45	5 mm	12 mm	44 mm



5mm Ripple Claw Forcep

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1226	Ratcheting Handle with Rotation	5 mm	23 mm	36 cm
98-1227	Slide Lock Handle	5 mm	23 mm	34 cm
98-1228	Slide Lock Handle with Rotation	5 mm	23 mm	36 cm



5mm Mini-Micro Dissecting Forcep

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1229	Monopolar Handle	5 mm	12 mm	35 cm
98-1230	Monopolar Handle with Rotation	5 mm	12 mm	35 cm
98-1231	Standard Handle with Rotation	5 mm	12 mm	35 cm

New Products

Symmetry Surgical Laparoscopic Instruments




5mm Flat Dissector

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
90-1040	Monopolar Handle	5 mm	13 mm	35 cm
90-1040L	Monopolar Handle, Opti-Length®	5 mm	13 mm	45 cm



5mm Biopsy Scoop Forcep

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
90-1075	Monopolar Handle	5 mm	13 mm	35 cm
90-1075L	Monopolar Handle, Opti-Length®	5 mm	13 mm	45 cm



5mm Tenaculum

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1232	Auto Closure Handle	5 mm	23 mm	34 cm
98-1232L	Auto Closure Handle, Opti-Length®	5 mm	23 mm	45 cm
98-1233	Auto Closure Handle with Rotation	5 mm	23 mm	36 cm
98-1234	Monopolar Handle with Rotation	5 mm	23 mm	36 cm
98-1235	Ratcheting Handle	5 mm	23 mm	36 cm
98-1236	Ratcheting Handle with Rotation	5 mm	23 mm	36 cm
98-1236L	Ratcheting Handle with Rotation, Opti-Length®	5 mm	23 mm	45 cm
98-1237	Slide Lock Handle	5 mm	23 mm	36 cm
98-1238	Slide Lock Handle with Rotation	5 mm	23 mm	36 cm
98-1238L	Slide Lock Handle with Rotation, Opti-Length®	5 mm	23 mm	45 cm
98-1239	Standard Handle with Rotation	5 mm	23 mm	36 cm

New Products

Symmetry Surgical Laparoscopic Instruments



5mm Medium Curved Dissector

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1240	Monopolar Handle with Rotation	5 mm	17 mm	36 cm
98-1240L	Monopolar Handle with Rotation, Opti-Length®	5 mm	17 mm	45 cm



Fenestrated Bowel Grasper, 5mm

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1241	Monopolar Handle with Rotation	5 mm	20 mm	36 cm
98-1241L	Monopolar Handle with Rotation, Opti-Length®	5 mm	20 mm	45 cm
98-1242	Ratcheting Handle	5 mm	20 mm	36 cm
98-1243	Ratcheting Handle with Rotation	5 mm	20 mm	36 cm
98-1244	Slide Lock Handle with Rotation	5 mm	20 mm	36 cm
98-1244L	Slide Lock handle with Rotation, Opti-Length®	5 mm	20 mm	45 cm
98-1245	Standard Handle with Rotation	5 mm	20 mm	36 cm
98-1245L	Standard Handle with Rotation, Opti-Length®	5 mm	20 mm	45 cm



5mm Curved DeBakey Tissue Forcep

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1246	Monopolar Handle with Rotation	5 mm	27 mm	37 cm
98-1247	Ratcheting Handle with Rotation	5 mm	27 mm	37 cm
98-1248	Slide Lock Handle with Rotation	5 mm	27 mm	37 cm
98-1249	Standard Handle with Rotation	5 mm	27 mm	37 cm

New Products

Symmetry Surgical Laparoscopic Instruments



5mm Andrew Babcock Forcep

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1250	Monopolar Handle with Rotation	5 mm	19 mm	36 cm
98-1251	Ratcheting Handle with Rotation	5 mm	19 mm	36 cm
94-1089SL	Slide Lock Handle	5 mm	19 mm	36 cm
98-1252	Slide Lock Handle with Rotation	5 mm	19 mm	36 cm
98-1252L	Slide Lock Handle with Rotation, Opti-Length®	5 mm	19 mm	45 cm




5mm Maryland with Cross Serrations

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1253	Monopolar Handle with Rotation	5 mm	19 mm	36 cm
98-1253L	Monopolar Handle with Rotation, Opti-Length®	5 mm	19 mm	45 cm
98-1254	Ratcheting Handle	5 mm	19 mm	36 cm



Murray Grasper, 5mm

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1255	Ratcheting Handle	5 mm	11 mm	35 cm
97-1078RC	Ratcheting Handle with Rotation	5 mm	11 mm	35 cm
98-1256	Slide Lock Handle with Rotation	5 mm	11 mm	35 cm



5mm Double Action Spoon Forcep

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1257	Standard Handle	5 mm	13 mm	36 cm
98-1258	Standard Handle with Rotation	5 mm	13 mm	35 cm
98-1258L	Standard Handle with Rotation, Opti-Length®	5 mm	13 mm	44 cm

New Products

Symmetry Surgical Laparoscopic Instruments



5mm Meah Grasper

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1259	Ratcheting Handle with Rotation	5 mm	21 mm	36 cm
98-1260	Slide Lock Handle with Rotation	5 mm	21 mm	36 cm
98-1260L	Slide Lock Handle with Rotation, Opti-Length®	5 mm	21 mm	45 cm



10mm Angled Circle Grasper

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1261	Ratcheting Handle with Rotation	10 mm	24 mm	36 cm
98-1261L	Ratcheting Handle with Rotation, Opti-Length®	10 mm	24 mm	45 cm



10mm Glassman

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1262	Ratcheting Handle with Rotation	10 mm	60 mm	37 cm
98-1262L	Ratcheting Handle with Rotation, Opti-Length®	10 mm	60 mm	45 cm



10mm Right Angle Dissector

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1263	Ratcheting Handle	10 mm	37 mm	37 cm
98-1264	Standard Handle	10 mm	37 mm	37 cm
98-1265	Standard Handle with Rotation	10 mm	37 mm	37 cm

New Products

Symmetry Surgical Laparoscopic Instruments



10mm Mixer Dissector

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1266	Standard Handle	10 mm	27 mm	37 cm
98-1267	Standard Handle with Rotation	10 mm	27 mm	37 cm



Claw Forcep, 10mm

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
90-1015	Auto Closure Handle	10 mm	29 mm	35 cm
98-1268L	Auto Closure Handle, Opti-Length®	10 mm	29 mm	45 cm
98-1269	Ratcheting Handle	10 mm	29 mm	35 cm
98-1270	Ratcheting Handle with Rotation	10 mm	29 mm	35 cm
98-1270L	Ratcheting Handle with Rotation, Opti-Length®	10 mm	29 mm	45 cm
98-1271	Standard Handle	10 mm	29 mm	35 cm
98-1272	Standard Handle with Rotation	10 mm	29 mm	35 cm



10mm Spoon Forcep

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1273	Auto Closure Handle	10 mm	25 mm	35 cm
98-1274	Auto Closure Handle with Rotation	10 mm	25 mm	35 cm
98-1273L	Auto Closure Handle, Opti-Length®	10 mm	25 mm	45 cm
90-1070	Ratcheting Handle	10 mm	25 mm	35 cm
90-1070L	Ratcheting Handle, Opti-Length®	10 mm	25 mm	45 cm
98-1275	Ratcheting Handle with Rotation	10 mm	25 mm	35 cm
98-1276	Standard Handle	10 mm	25 mm	35 cm
98-1277	Standard Handle with Rotation	10 mm	25 mm	35 cm
98-1277L	Standard Handle with Rotation, Opti-Length®	10 mm	25 mm	45 cm

New Products

Symmetry Surgical Laparoscopic Instruments



10mm Babcock, Short Jaw

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1278	Auto Closure with Rotation	10 mm	25 mm	35 cm
97-1083RC	Ratcheting Handle with Rotation	10 mm	25 mm	35 cm



10mm Babcock with Straight Serrations

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
97-1097RC	Ratchet Handle with Rotation	10 mm	35 mm	37 cm
98-1288L	Ratchet Handle with Rotation, Opti-Length®	10 mm	35 mm	45 cm



10mm Paddle Babcock

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1279	Ratchet Handle with Rotation	10 mm	35 mm	37 cm
98-1280L	Ratchet Handle with Rotation, Opti-Length®	10 mm	35 mm	45 cm



10mm Tenaculum

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
98-1281	Auto Closure Handle	10 mm	36 mm	38 cm
98-1282	Auto Closure Handle with Rotation	10 mm	36 mm	38 cm
98-1283	Ratcheting Handle	10 mm	36 mm	38 cm
98-1284	Ratcheting Handle with Rotation	10 mm	36 mm	38 cm
98-1284L	Ratcheting Handle with Rotation, Opti-Length®	10 mm	36 mm	45 cm
98-1285	Standard Handle with Rotation	10 mm	36 mm	38 cm

Accessories & Miscellaneous



Fan Retractor

CODE	DESCRIPTION	DIAMETER	JAW SIZE	LENGTH
24-1800	Pressure Lock Handle	10 mm	60 mm	24 cm
32-1800	Pressure Lock Handle	10 mm	60 mm	32 cm
38-1800	Pressure Lock Handle	10 mm	60 mm	38 cm
45-1800	Pressure Lock Handle	10 mm	60 mm	45 cm
24-1800R	Ratchet Lock Handle	10 mm	60 mm	24 cm
32-1800R	Ratchet Lock Handle	10 mm	60 mm	32 cm
38-1800R	Ratchet Lock Handle	10 mm	60 mm	38 cm
45-1800R	Ratchet Lock Handle	10 mm	60 mm	45 cm



Bariatric Hook


CODE	DESCRIPTION	DIAMETER	LENGTH
2010101	In-Line Handle	5 mm	45 cm



Bariatric Knot Pusher

CODE	DESCRIPTION	DIAMETER	LENGTH
2010201	In-Line Handle	10 mm	45 cm

Accessories & Miscellaneous



504805

Bookwalter® Universal Flex Arm System

Bookwalter® Table Post and Bookwalter® Coupler Required - Not Included

CODE	DESCRIPTION
50-4805	Universal Flex Arm System (Includes Flex Arm and 1 Each of the 4 Adapters)
50-4806	Ring Adapter for Universal Flex Arm
50-4807	Lap Instrument Adapter for Universal Flex Arm
50-4808	Camera Lens Adapter for Universal Flex Arm
50-4809	Camera Lens Adapter for Universal Flex Arm

Bookwalter® Laparoscopic Support Set

CODE	DESCRIPTION	QTY
50-5730	Laparoscopic Support Set Includes the Following:	
50-4683	Bookwalter® Table Post, 8 in	1
50-4554	Bookwalter® Post Coupling	1
50-4805	Bookwalter® Universal Flex Arm, Includes 4 Adapters Listed Below	1
50-4806	Bookwalter® Flex Arm Ring Adapter	1
50-4607	Bookwalter® Flex Arm Laparoscopic Instrument Adapter	1
50-4608	Bookwalter® Flex Arm Camera Lens Adapter	1
50-4609	Bookwalter® Flex Arm Camera Housing Adapter	1
68-1593	Fan Retractor, 5mm, Articulating, 3 Fingers	1
68-1594	Fan Retractor, 10mm, Articulating, 5 Fingers	1



Nathanson Flex Arm

Holds 5 mm and 10 mm Instruments, For Use with Bookwalter® Table Post and Bookwalter® Coupler

CODE	DESCRIPTION
SPO1115-45	Flex Arm

Accessories & Miscellaneous

Circular Anvil Grasper

CODE	DESCRIPTION	DIAMETER	LENGTH
2010202	In-Line Handle	10 mm	45 cm



Tip

38-5411

Dodd Knot Pusher

CODE	DESCRIPTION	LENGTH
38-5401	0-Suture	38 cm
38-5411	2-0 Suture	38 cm



Fan Retractors

To Be Used with the Universal Flex Arm on LAP-23

CODE	DESCRIPTION	DIAMETER
681590	Fan Liver Retractor, Three Fingers	5 mm
681593	Fan Liver Retractor, Articulating, Three Fingers	5 mm
681591	Fan Liver Retractor, Articulating, Five Fingers	10 mm
681594	Fan Liver Retractor, Articulating, Five Fingers	10 mm



Tip

Gastric Band Passer

CODE	DESCRIPTION	DIAMETER	LENGTH
2010301	In-Line Handle	15 mm	45 cm

Accessories & Miscellaneous



32-5397

Handle

To Be Used with 12 Series Reposable Scissors

CODE	DESCRIPTION	DIAMETER
32-5399	Reposable, Monopolar, Non-Locking with Metzenbaum Tip	32 cm
32-5397	Reposable, Monopolar, Non-Locking	32 cm
38-5397	Reposable, Monopolar, Non-Locking	38 cm



90-6004

Hasson "S" Retractor

Set of 2

CODE	WIDTH		DIAMETER	LENGTH
90-6003	¼ in	6.35 mm	6 in	152.4 mm
90-6004	½ in	12.7 mm	6 in	152.4 mm
90-6005	¾ in	10 mm	7 ½ in	180 mm
90-6006	5/8 in	15 mm	7 ½ in	180 mm
90-6007	¾ in	19 mm	7 ½ in	180 mm


TIP
12-5383

Hook Reposable Scissors

For Use with Reposable Handles 32-5397 and 38-5397

CODE	DESCRIPTION	TIP SIZE	DIAMETER	LENGTH
12-5383	Monopolar, Non-Locking	8 mm	5 mm	32 cm
12-5393	Monopolar, Non-Locking	8 mm	5 mm	38 cm

Accessories & Miscellaneous



Irrigating/Aspirating Cannula

CODE	DESCRIPTION
90-5070	Dual Valve



Luer Lock Cap

CODE	DESCRIPTION
90-0000	Choli, 10/Pack



TIP
12-5382

Metzenbaum Repposable Scissors

For Use with Repposable Handles 32-5397, 32-5399 and 38-5397

CODE	DESCRIPTION	DIAMETER	TIP SIZE	QTY
12-5382	Curved, 45° Opening	5 mm	11 mm	5/Pack
12-5392	Curved, 45° Opening	5 mm	11 mm	
12-5384	Curved, Delicate, 45° Opening	5 mm	11 mm	5/Pack
12-5394	Curved, Delicate, 45° Opening	5 mm	11 mm	

Accessories & Miscellaneous



Myoma Screw

CODE	DESCRIPTION	DIAMETER	LENGTH
90-5600	Screw	5 mm	32 cm
90-5605	Screw	5 mm	45 cm

Nathanson Retractor Hooks

CODE	DESCRIPTION	DIAMETER	DIMENSIONS
52420-00	Kit includes One Small, One Medium, One Large	6.5 mm	
52420-01	Small	6.5 mm	3 1/2 x 7 x 4 1/2 in 89 x 178 x 114 mm
52420-02	Medium	6.5 mm	4 1/2 x 8 1/2 x 8 in 114 x 216 x 203 mm
52420-03	Large	6.5 mm	6 1/2 x 11 x 10 in 165 x 279 x 254 mm

90-5600



Reducer Sleeve

CODE	SIZE
90-5030	5-4 mm

TIP
12-5391

Reusable Scissors

Use with Reusable Handle 32-5399


CODE	DESCRIPTION	DIAMETER	TIP SIZE
12-5391	Curved, Single, 45° Opening	5 mm	8 mm

Accessories & Miscellaneous



Roeder Loop Applicator


CODE	DESCRIPTION	DIAMETER	LENGTH
91-0630	Loop Applicator	5 mm	15 cm



Tip
29-5403

Saye/Reddick Knot Pusher

CODE	DESCRIPTION	LENGTH
29-5403	Closed	29 cm
38-5403	Closed	38 cm
38-5402	Opened	38 cm



Tip


38-5506

Saye/Reddick Needle Holder

Pressure Lock Handle, with Flush Port

CODE	DESCRIPTION	DIAMETER	JAW SIZE	LENGTH
38-5505	Straight, Double-Action	5 mm	11 mm	38 cm
38-5506	Straight, Single-Action	5 mm	11 mm	38 cm
45-5506	Straight, Single-Action	5 mm	11 mm	45 cm
24-5507	Curved, Single-Action	5 mm	11 mm	24 cm
32-5507	Curved, Single-Action	5 mm	11 mm	32 cm
38-5507	Curved, Single-Action	5 mm	11 mm	38 cm

Accessories & Miscellaneous



Trocars

Hasson Open Cannula

CODE	DESCRIPTION	DIAMETER
90-9101	With Suture Tying Posts	11 mm
90-9102	With Suture Tying Posts	10 mm



Ultra™ Ergonomic Needle Holder


CODE	DESCRIPTION	DIAMETER	LENGTH
90-3035	Ring Handle In-Line with Shaft, Ratchet Handle, Single-Action	5 mm	32 cm



Ultra™ Knot Pusher

CODE	DESCRIPTION	DIAMETER	LENGTH
90-5099	Opened	5 mm	32 cm

Accessories & Miscellaneous



Ultra™ Needle Holder

Axial Handle

CODE	DESCRIPTION	DIAMETER	LENGTH
90-2011	Needle Holder, Curved Left, Tungsten Carbide	5 mm	35 cm



90-5053

Ultra™ Palpation Probe

CODE	DESCRIPTION	DIAMETER	LENGTH
90-5053	With Centimeter Graduations	5 mm	32 cm
90-5053L	With Centimeter Graduations	5 mm	45 cm



90-5049

Ultra™ Puncture Needle with Aspiration

Puncture Needle Aspirator

CODE	DESCRIPTION	DIAMETER	LENGTH
90-5048	10 Gauge	5 mm	32 cm
90-5049	12 Gauge	5 mm	32 cm
90-5051	15 Gauge	5 mm	32 cm
90-5050	19 Gauge	5 mm	32 cm
90-5052	22 Gauge	5 mm	32 cm
90-5051L	15 Gauge	5 mm	45 cm
90-5050L	19 Gauge	5 mm	45 cm
90-5052L	22 Gauge	5 mm	45 cm

Accessories & Miscellaneous



Ultra™ Tungsten Carbide Needle Holder

90-3038

Flush Port, Squeeze Handle with Ratchet

CODE	DESCRIPTION	DIAMETER	LENGTH
90-3039	Straight	5 mm	32 cm
90-3038	Curved	5 mm	32 cm
90-3038L	Curved	5 mm	45 cm



Verres Needle

90-5054

CODE	DESCRIPTION	DIAMETER	LENGTH
90-5058	Needle	2 mm	7.5 cm
90-5054	Needle	2 mm	10 cm
90-5055	Needle	2 mm	12.5 cm
90-5059	Needle	2 mm	20 cm

Cables



90-1020

Reusable Bipolar Cable

CODE	DESCRIPTION	LENGTH
90-1020	For Symmetry/Wolf Kleppinger Type Forceps and Wolf Generator	10 ft
90-1021	For Symmetry/Wolf Kleppinger Forceps and Birtcher, Cabot and Valley Lab Type Generator	10 ft

Dissectors

Secto® Laparoscopic Dissectors

5 mm x 400 mm Stick, 10 per box



CODE	DESCRIPTION
81-1003	Secto® Cherry Laparoscopic Dissector, 10 mm Cherry Sponge Spun Cotton Tip
81-1004	Secto® Bullet Laparoscopic Dissector, 5 mm Bullet Sponge Spun Cotton Tip
81-1005	Secto® Kittner Laparoscopic Dissector, ¼ in x ⅝ in Kittner Sponge Dissector Tip, Hand-Sewn Cotton Chevron Tape Roll


MultiPak

MultiPak Accessories

ITEM	CODE	DESCRIPTION	SIZE
	5200	Endoscopic Silicone Bar Holder, for Wire Rack 5100-5113, 7113	7 in 177 mm
	5200N	Endoscopic Silicone Bar Holder, for Wire Rack 5108 or 5109	6.5 in 165 mm

*Silicone accessories include plastic holding bar and screws; excluding 5310


ScopePak



ScopePak Cases

CODE	DESCRIPTION	DIMENSIONS (INTERNAL)
9417-1	OnePak, 17 in, Configured to Accommodate 1 Scope	3 ⅛ x 1 7/10 x 1 ½ in 82.5 x 43 x 38 mm


ScopePak



ScopePak Base & Lid Only

9321-00

CODE	DESCRIPTION	DIMENSIONS (INTERNAL)	
9312-00	Small Base and Lid	3 x 11 $\frac{1}{10}$ x 1 $\frac{3}{10}$ in	76 x 282 x 33 mm
9317-00	Medium Base and Lid	3 x 16 x 1 $\frac{3}{10}$ in	76 x 406 x 33 mm
9321-00	Large Base and Lid	3 x 20 x 1 $\frac{3}{10}$ in	76 x 508 x 33 mm
9326-00	Extra-Large Base and Lid	4 x 26 x 2 $\frac{4}{5}$ in	102 x 660 x 71 mm
9331-00	Long Base and Lid	4 $\frac{7}{10}$ x 30 $\frac{3}{10}$ x 2 $\frac{4}{5}$ in	119 x 769 x 71 mm



ScopePak Base & Lid with Scope Holder Sets

9312


CODE	DESCRIPTION	INCLUDES	DIMENSIONS (INTERNAL)	
9312	Small Base and Lid	2/Each 5350 1/Each 5351	3 x 11 $\frac{1}{10}$ x 1 $\frac{3}{10}$ in	76 x 282 x 33 mm
9317	Medium Base and Lid	4/Each 5353 2/Each 5352	3 x 16 x 1 $\frac{3}{10}$ in	76 x 406 x 33 mm
9321	Large Base and Lid	2/Each 5350 1/Each 5351	3 x 20 x 1 $\frac{3}{10}$ in	76 x 508 x 33 mm
9326	Extra-Large Base and Lid	1/Each 5350 1/Each 5351	4 x 26 x 2 $\frac{4}{5}$ in	102 x 660 x 71 mm
9331	Long Base and Lid	4/Each 5353 2/Each 5352	4 $\frac{7}{10}$ x 30 $\frac{3}{10}$ x 2 $\frac{4}{5}$ in	119 x 769 x 71 mm

Symmetry Surgical Laparoscopic Instruments

FEATURE: PUSH DRIVE ACTUATION

Patented link mechanism eliminates the need for wings at the distal tip, while allowing for strong tactile feedback.

- Flushport to Assist in Cleaning
- Internal One-Way Seal helps prevent bio burden from flowing back up the inner lumens
- Ergonomic Handle Design for comfortable fit with precise control
- Laparoscopic Instruments Available In 45 cm Length
- Made in the USA



Symmetry Surgical Laparoscopic Instruments



Aggressive Retention Grasping Forceps
60° Opening

32-5101


CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
32-5101	Pressure Lock Monopolar Handle	5 mm	15 mm	32 cm
38-5101	Pressure Lock Monopolar Handle	5 mm	15 mm	38 cm
32-5101R	Ratchet Lock Handle	5 mm	15 mm	32 cm
38-5101R	Ratchet Lock Handle	5 mm	15 mm	38 cm



Allis Grasping Forceps
60° Opening

32-5103

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
32-5103	Pressure Lock Handle	5 mm	18 mm	32 cm
38-5103	Pressure Lock Handle	5 mm	18 mm	38 cm
38-5103R	Ratchet Lock Handle	5 mm	18 mm	38 cm



Allis Grasping Forceps with DeBakey Teeth
60° Opening

32-5110

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
24-5110	Pressure Lock Handle	5 mm	18 mm	24 cm
32-5110	Pressure Lock Handle	5 mm	18 mm	32 cm

Symmetry Surgical Laparoscopic Instruments



32-5117

Atraumatic Fenestrated Grasping Forceps

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
32-5117	Pressure Lock, Monopolar Handle, 60° Opening	5 mm	20 mm	32 cm
38-5117	Pressure Lock, Monopolar Handle, 60° Opening	5 mm	20 mm	38 cm



24-5157

Atraumatic Fenestrated Parallel Jaw Grasping Forceps

For Enhanced Grasping Power without the Need for Excessive Closing Force

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
24-5157	Pressure Lock Handle	5 mm	16 mm	24 cm
32-5157	Pressure Lock Handle	5 mm	16 mm	32 cm
38-5157	Pressure Lock Handle	5 mm	16 mm	38 cm
45-5157	Pressure Lock Handle	5 mm	16 mm	45 cm
24-5158	Pressure Lock Handle	5 mm	25 mm	24 cm
32-5158	Pressure Lock Handle	5 mm	25 mm	32 cm
38-5158	Pressure Lock Handle	5 mm	25 mm	38 cm
45-5158	Pressure Lock Handle	5 mm	25 mm	45 cm
24-5157R	Ratchet Lock Handle	5 mm	16 mm	24 cm
38-5157R	Ratchet Lock Handle	5 mm	16 mm	38 cm
45-5157R	Ratchet Lock Handle	5 mm	16 mm	45 cm
24-5158R	Ratchet Lock Handle	5 mm	25 mm	24 cm
32-5158R	Ratchet Lock Handle	5 mm	25 mm	32 cm
38-5158R	Ratchet Lock Handle	5 mm	25 mm	38 cm
45-5158R	Ratchet Lock Handle	5 mm	25 mm	45 cm

Symmetry Surgical Laparoscopic Instruments



Atraumatic Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
38-1600R	Grasper Locking Handle	10 mm	38 cm



Babcock Grasping Forceps

60° Opening

24-5104

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
24-5104	Pressure Lock Handle	5 mm	18 mm	24 cm
32-5104	Pressure Lock Handle	5 mm	18 mm	32 cm
38-5104	Pressure Lock Handle	5 mm	18 mm	38 cm
45-5104	Pressure Lock Handle	5 mm	18 mm	45 cm



32-5109

Babcock Grasping Forceps with DeBaKey Teeth

60° Opening

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
32-5109	Pressure Lock Handle	5 mm	18 mm	32 cm
32-1003	Pressure Lock Handle	5 mm	23 mm	32 cm
32-5109R	Ratchet Lock Handle	5 mm	18 mm	32 cm
32-1003R	Ratchet Lock Handle	5 mm	23 mm	38 cm
38-5109R	Ratchet Lock Handle	5 mm	18 mm	38 cm
38-1003R	Ratchet Lock Handle	5 mm	23 mm	38 cm

Symmetry Surgical Laparoscopic Instruments



DeBakey Grasping Forceps

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
32-5610	Pressure Lock Handle	5 mm	10 mm	32 cm
45-5610	Pressure Lock Handle	5 mm	10 mm	45 cm
24-5620	Pressure Lock Handle	5 mm	20 mm	24 cm
32-5620	Pressure Lock Handle	5 mm	20 mm	32 cm
24-5640	Pressure Lock Handle	5 mm	40 mm	24 cm
32-5640	Pressure Lock Handle	5 mm	40 mm	32 cm
38-5640	Pressure Lock Handle	5 mm	40 mm	38 cm
45-5640	Pressure Lock Handle	5 mm	40 mm	45 cm
38-5640R	Ratchet Lock Handle	5 mm	10 mm	38 cm
45-5610R	Ratchet Lock Handle	5 mm	10 mm	45 cm
38-5620R	Ratchet Lock Handle	5 mm	20 mm	38 cm



Dissecting Forceps, Left

Non-Locking, Monopolar Handle, 90° Opening

CODE	DIAMETER	TIP SIZE	LENGTH
245201	5 mm	11 mm	24 cm
32-5201	5 mm	11 mm	32 cm
38-5201	5 mm	11 mm	38 cm



Dissecting Forceps, Soft Left

Non-Locking, Monopolar Handle, 90° Opening

CODE	DIAMETER	TIP SIZE	LENGTH
24-5202A	5 mm	18 mm	24 cm
32-5202	5 mm	18 mm	32 cm
38-5202	5 mm	18 mm	38 cm

Symmetry Surgical Laparoscopic Instruments



38-5205

Dissecting Forceps, Snub Nose

Non-Locking, Monopolar Handle, 90° Opening

CODE	DIAMETER	TIP SIZE	LENGTH
24-5205	5 mm	11 mm	24 cm
32-5205	5 mm	11 mm	32 cm
38-5205	5 mm	11 mm	38 cm



24-5200

Dissecting Forceps, Straight

Non-Locking, Monopolar Handle

CODE	DIAMETER	TIP SIZE	LENGTH
24-5200	5 mm	13 mm	24 cm
32-5200	5 mm	13 mm	32 cm
38-5200	5 mm	13 mm	38 cm



Tip
32-5203

Dissecting Forceps, Tapered

Non-Locking, Monopolar Handle, 90° Opening

CODE	DIAMETER	TIP SIZE	LENGTH
24-5203A	5 mm	15 mm	24 cm
32-5203	5 mm	15 mm	32 cm

Symmetry Surgical Laparoscopic Instruments



Duckbill Dissecting Forceps

Monopolar Handle, 90° Opening

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
32-5204	Non-Locking Handle	5 mm	11 mm	32 cm



Duval Lung Forceps

With DeBakey Teeth, 60° Opening

24-1500

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
24-1500	Pressure Lock Handle	10 mm	53 mm	24 cm
32-1500	Pressure Lock Handle	10 mm	53 mm	32 cm


Gemini Dissecting Forceps

Monopolar Handle, 90° Opening



CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
32-1206	Non-Locking Handle	10 mm	18 mm	32 cm

Symmetry Surgical Laparoscopic Instruments



Glassman Grasping Forceps


60° Opening

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
32-5720	Pressure Lock Handle	5 mm	20 mm	32 cm



Grasping Forceps

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
32-5207	Non-Locking, Monopolar Handle, Tapered	5 mm	15 mm	32 cm



32-5108R

Grasping Forceps

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
32-5108	Pressure Lock, Monopolar Handle, Short Nose, 60° Opening	5 mm	11 mm	32 cm
32-5108R	Ratchet Lock Handle, Short Nose, 60° Opening	5 mm	11 mm	32 cm
38-5108R	Ratchet Lock Handle, Short Nose, 60° Opening	5 mm	11 mm	38 cm

Symmetry Surgical Laparoscopic Instruments



32-5107R

Grasping Forceps

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
32-5107	Pressure Lock, Monopolar Handle, Fenestrated	5 mm	18 mm	32 cm
38-5107	Pressure Lock, Monopolar Handle, Fenestrated	5 mm	18 mm	38 cm
32-5107R	Ratchet Handle, Fenestrated	5 mm	18 mm	32 cm



32-1001

Grasping & Extracting Forceps

2x3 Teeth

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
32-1001	Pressure Lock Handle	10 mm	27 mm	32 cm
38-1001	Pressure Lock Handle	10 mm	27 mm	38 cm
32-1001R	Ratchet Handle	10 mm	27 mm	32 cm
32-5001R	Ratchet Lock Handle	5 mm	27 mm	32 cm
38-5001R	Ratchet Lock Handle	5 mm	27 mm	38 cm



32-5303

Hook Scissors

Single-Action, 45° Opening

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
32-5303	Non-Locking, Monopolar Handle	5 mm	9 mm	32 cm
38-5303	Non-Locking, Monopolar Handle	5 mm	9 mm	38 cm
45-5303	Non-Locking, Monopolar Handle	5 mm	9 mm	45 cm

Symmetry Surgical Laparoscopic Instruments



24-5302

Metzenbaum Scissors, Curved

Monopolar Handle, 45° Opening

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
24-5302	Non-Locking Handle	5 mm	11 mm	24 cm
32-5302	Non-Locking Handle	5 mm	11 mm	32 cm
38-5302	Non-Locking Handle	5 mm	11 mm	38 cm
45-5302	Non-Locking Handle	5 mm	11 mm	45 cm
32-5304	Non-Locking Handle, Delicate	5 mm	11 mm	32 cm
38-5304	Non-Locking Handle, Delicate	5 mm	11 mm	38 cm



Retention Grasping Aggressive Forceps

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
32-5106	Pressure Lock, Monopolar Handle, Tapered, 60° Opening	5 mm	15 mm	32 cm



24-5100

Retention Grasping Forceps

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
24-5100	Pressure Lock, Monopolar Handle	5 mm	15 mm	24 cm
27-5100	Pressure Lock, Monopolar Handle	5 mm	15 mm	27 cm
32-5100	Pressure Lock, Monopolar Handle	5 mm	15 mm	32 cm
32-5100R	Ratchet Lock Handle	5 mm	15 mm	32 cm
38-5100R	Ratchet Lock Handle	5 mm	15 mm	38 cm

Symmetry Surgical Laparoscopic Instruments



Scissors

Single-Action

38-5300

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
32-5300	Straight, Non-Locking, Monopolar Handle	5 mm	8 mm	32 cm
38-5300	Straight, Non-Locking, Monopolar Handle	5 mm	8 mm	38 cm
24-5301	Curved, Non-Locking, Monopolar Handle	5 mm	8 mm	24 cm
32-5301	Curved, Non-Locking, Monopolar Handle	5 mm	8 mm	32 cm



Snake Tooth (COBRA) Grasping Forceps

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
32-5102	Pressure Lock, Monopolar Handle, 60° Opening	5 mm	15 mm	32 cm



Stone Retrieval Forceps

32-5900

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
32-5900	Non-Locking Handle, Non-Cautery	5 mm	11 mm	32 cm
32-1900	Non-Locking Handle, Non-Cautery	10 mm	27 mm	32 cm
38-5900	Non-Locking Handle, Non-Cautery	5 mm	11 mm	38 cm

Symmetry Surgical Laparoscopic Instruments



Wave Atraumatic Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
38-5209	Pressure Lock Handle	5 mm	38 cm
45-5209	Pressure Lock Handle	5 mm	45 cm



Wave Atraumatic Grasping Forceps

Single-Action, 45° Opening

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
24-5127	Pressure Lock, Monopolar Handle	5 mm	20 mm	24 cm
32-5127	Pressure Lock, Monopolar Handle	5 mm	20 mm	32 cm
38-5127	Pressure Lock, Monopolar Handle	5 mm	20 mm	38 cm
32-5137	Pressure Lock Handle, Long, Square Teeth	5 mm	20 mm	32 cm
32-5147	Pressure Lock Handle, Long, Square Teeth	5 mm	29 mm	32 cm
38-5127	Pressure Lock Handle	5 mm	20 mm	38 cm
38-5147	Pressure Lock Handle	5 mm	29 mm	38 cm
32-5127R	Ratchet Lock Handle	5 mm	20 mm	32 cm
32-5137R	Ratchet Lock Handle	5 mm	20 mm	32 cm
32-5147R	Ratchet Lock Handle	5 mm	29 mm	32 cm
38-5147R	Ratchet Lock Handle	5 mm	29 mm	38 cm

Symmetry Surgical Laparoscopic Instruments




Allis Tissue Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1077	Ratchet Handle	10 mm	32 cm



Babcock Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1079	Ratchet Handle, Double-Action	10 mm	32 cm



Biopsy Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1075	Monopolar, Fully-Insulated, Single-Action	5 mm	32 cm
90-1075L	Monopolar, Fully-Insulated, Single-Action	5 mm	45 cm

Symmetry Surgical Laparoscopic Instruments



Bottle Nose (Dolphin Nose) Dissecting Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1030	Monopolar, Fully-Insulated, Double-Action	5 mm	32 cm



Cholangiogram Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
91-5035	Ratchet Handle, Accepts Catheter Up to 7 French, Double-Action	5 mm	39 cm
90-1076	Cholangiogram Guide, Curved		32 cm



Claw Grasping Extractor Forceps

2 x 3 Heavy Teeth, Single-Action

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1005	Spring Handle	5 mm	33 cm
90-1015	Spring Handle	10 mm	33 cm
90-1015L	Spring Handle	10 mm	45 cm

Symmetry Surgical Laparoscopic Instruments



Cup Tooth Dissecting Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1035	Non-Locking, Monopolar, Fully-Insulated, Double-Action	5 mm	32 cm
90-1035L	Non-Locking, Monopolar, Fully-Insulated, Double-Action	5 mm	45 cm



Cup Tooth Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1055	Ratchet Handle, Double-Action	5 mm	32 cm



Duckbill Dissecting Forceps

Flattened Serrated Jaw

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1040	Monopolar, Fully Insulated	5 mm	32 cm

Symmetry Surgical Laparoscopic Instruments



Duval Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1083	Ratchet Handle	10 mm	32 cm



Fundus Trimax Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1050	Ratchet Handle, Aggressive Jaws, Double-Action	5 mm	32 cm



Hook Biopsy Forceps

CODE	DIAMETER	LENGTH
90-1073	5 mm	33 cm

Symmetry Surgical Laparoscopic Instruments



Hook Scissors

Non-Locking

CODE	DESCRIPTION	DIAMETER	LENGTH
90-2020	Monopolar, Fully-Insulated, Single-Action	5 mm	32 cm
90-2020L	Monopolar, Fully-Insulated, Single-Action	5 mm	45 cm



Maryland Dissecting Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1067	Monopolar, Fully-Insulated, Double-Action	5 mm	32 cm
90-1065	Monopolar, Fully-Insulated, Double-Action	10 mm	32 cm



Micro Fenestrated Dissecting Forceps

Fine Serrations

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1044	Monopolar, Fully-Insulated	5 mm	32 cm

Symmetry Surgical Laparoscopic Instruments



Mixer Dissecting Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1068	Non-Locking Handle, Monopolar, Fully Insulated	10 mm	32 cm



Pennington Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1082	Ratchet Handle	10 mm	32 cm



Retrograde Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1031	Ratchet Handle	5 mm	32 cm

Symmetry Surgical Laparoscopic Instruments



Grasper

CODE	DESCRIPTION	DIAMETER	LENGTH
90-2210	Insulated, Monopolar, Single-Action, Fine Tungsten Carbide Lower Blade, Serrated	5 mm	33 cm



Snake Tooth (COBRA) Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1069	Ratchet Handle, 2 x 3 Teeth, Double-Action	5 mm	32 cm



Stone Retrieval Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1070	Spring Handle, Single-Action	10 mm	32 cm

Symmetry Surgical Laparoscopic Instruments




Tenaculum Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1008	Single Tooth, Double-Action	5 mm	32 cm



Tooth Cup Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1072	Monopolar, Non-Locking	5 mm	32 cm



Tungsten Carbide Needle Holder

CODE	DESCRIPTION	DIAMETER	LENGTH
90-2010	Straight In-Line Handle	5 mm	35 cm

Symmetry Surgical Laparoscopic Instruments

360° rotation with flush port and Direct Drive system.
Laparoscopic instrumentation available in 45 cm length (-45).



94-1078L

Allis Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
94-1078L	Direct Drive, Locking Ratchet on Demand	5 mm	33 cm
94-1078L-45	Direct Drive, Locking Ratchet on Demand	5 mm	45 cm
94-1088L	Direct Drive, Locking Ratchet on Demand	5 mm	33 cm
94-1078SL	Direct Drive, Slide Lock Handle	5 mm	33 cm
94-1078SL-45	Direct Drive, Slide Lock Handle	5 mm	45 cm




94-1053L

Atraumatic Grasping Forceps

Monopolar Handle

CODE	DESCRIPTION	DIAMETER	LENGTH
94-1053L	Direct Drive, Locking Ratchet on Demand, Double-Action	5 mm	33 cm
94-1053SL	Direct Drive, Slide Lock Handle, Double-Action	5 mm	33 cm
94-1053SL-45	Slide Lock Handle, Double-Action	5 mm	45 cm



94-1089L

Babcock Grasping Forceps with DeBakey Teeth

CODE	DESCRIPTION	DIAMETER	LENGTH
94-1089L	Direct Drive, Locking Ratchet on Demand	5 mm	33 cm
94-1089SL	Direct Drive, Locking Ratchet on Demand	5 mm	33 cm

Symmetry Surgical Laparoscopic Instruments

360° rotation with flush port and Direct Drive system.
Laparoscopic instrumentation available in 45 cm length (-45).



Biopsy Punch

CODE	DESCRIPTION	DIAMETER	LENGTH
94-1075	Direct Drive, Non-Locking, Monopolar Handle	5 mm	33 cm



Bottle Nose (Dolphin Nose) Dissecting Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
94-1030	Direct Drive, Non-Locking, Double-Action	5 mm	33 cm



Bottle Nose (Dolphin Nose) Dissecting Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
94-1034	Direct Drive, Non-Locking, Monopolar Handle	5 mm	33 cm

Symmetry Surgical Laparoscopic Instruments


360° rotation with flush port and Direct Drive system.
Laparoscopic instrumentation available in 45 cm length (-45).



Claw Grasping Extractor Forceps

2x3 Heavy Teeth

CODE	DESCRIPTION	DIAMETER	LENGTH
94-1005L	Direct Drive, Locking Ratchet on Demand	5 mm	33 cm
94-1015L	Direct Drive, Locking Ratchet on Demand	10 mm	33 cm
94-1005L-45	Direct Drive, Locking Ratchet on Demand	5 mm	45 cm
94-1015L-45	Direct Drive, Locking Ratchet on Demand	10 mm	45 cm



Cup Tooth Dissecting Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
94-1035	Direct Drive, Non-Locking	5 mm	33 cm



Cup Tooth Grasping Forceps

Non-Rotating

CODE	DESCRIPTION	DIAMETER	JAW SIZE	LENGTH
94-1035SL	Direct Drive, Slide Lock Handle, Double-Action	5 mm	13 mm	33 cm
94-1035L	Direct Drive, Locking Ratchet on Demand	5 mm	13 mm	33 cm

Symmetry Surgical Laparoscopic Instruments

360° rotation with flush port and Direct Drive system.
Laparoscopic instrumentation available in 45 cm length (-45).



DeBakey Grasping Forceps

Non-Rotating

CODE	DESCRIPTION	DIAMETER	JAW SIZE	LENGTH
94-1091L	Direct Drive, Locking Ratchet on Demand	5 mm	20 mm	33 cm
94-1091L-45	Direct Drive, Locking Ratchet on Demand	5 mm	20 mm	45 cm
94-1091SL	Direct Drive, Slide Lock Handle	5 mm	20 mm	33 cm
94-1092SL	Direct Drive, Slide Lock Handle	5 mm	40 mm	33 cm



Dissecting Forceps

Soft Curve

CODE	DESCRIPTION	DIAMETER	LENGTH
94-1064	Direct Drive, Non-Locking, Monopolar Handle	5 mm	33 cm



Dissecting Forceps

Blunt Tip

CODE	DESCRIPTION	DIAMETER	LENGTH
94-1064	Direct Drive, Non-Locking, Monopolar Handle	5 mm	33 cm

Symmetry Surgical Laparoscopic Instruments

360° rotation with flush port and Direct Drive system.
Laparoscopic instrumentation available in 45 cm length (-45).



94-1032L

Tip

Grasping Forceps

Blunt Tip

CODE	DESCRIPTION	DIAMETER	LENGTH
94-1032L	Direct Drive, Locking Ratchet on Demand	5 mm	33 cm
94-1032L-45	Direct Drive, Locking Ratchet on Demand	5 mm	45 cm



Tip

Hook Scissors

CODE	DESCRIPTION	DIAMETER	LENGTH
94-2020	Direct Drive, Non-Locking	5 mm	33 cm



94-1093

Hunter Bowel Grasping Forceps

CODE	DESCRIPTION	DIAMETER	JAW SIZE	LENGTH
94-1093	Direct Drive, Non-Locking, Monopolar Handle	5 mm	20 mm	33 cm
94-1093-45	Direct Drive, Non-Locking, Monopolar Handle	5 mm	20 mm	45 cm
94-1093L	Direct Drive, Locking, Ratchet on Demand, Monopolar Handle	5 mm	20 mm	33 cm
94-1093L-45	Direct Drive, Locking, Ratchet on Demand, Monopolar Handle	5 mm	20 mm	45 cm

Symmetry Surgical Laparoscopic Instruments


360° rotation with flush port and Direct Drive system.
Laparoscopic instrumentation available in 45 cm length (-45).



Maryland Dissecting Forceps

Curved

CODE	DESCRIPTION	DIAMETER	LENGTH
94-1067-45	Direct Drive, Non-Locking, Monopolar Handle	5 mm	45 cm



Mixer Dissecting Forceps

Right Angle

CODE	DESCRIPTION	DIAMETER	LENGTH
94-1066	Direct Drive, Non-Locking, Monopolar Handles	5 mm	33 cm
94-1066-45	Direct Drive, Non-Locking, Monopolar Handles	5 mm	45 cm



Mixer Dissecting Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
94-1068	Direct Drive, Non-Locking	10 mm	33 cm

Symmetry Surgical Laparoscopic Instruments

360° rotation with flush port and Direct Drive system.
Laparoscopic instrumentation available in 45 cm length (-45).



Micro Fenestrated Dissecting Forceps

Fine Serration

CODE	DESCRIPTION	DIAMETER	LENGTH
94-1044	Direct Drive, Monopolar, Fully-Insulated	5 mm	33 cm



Needle Nose Dissecting Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
94-1045	Direct Drive, Non-Locking, Monopolar Handle	5 mm	33 cm



Retraction Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
94-1056L	Direct Drive, Locking Ratchet on Demand, Single-Action	5 mm	33 cm

Symmetry Surgical Laparoscopic Instruments


360° rotation with flush port and Direct Drive system.
Laparoscopic instrumentation available in 45 cm length (-45).



Slide Lock Retraction Grasping Forceps

Non-Rotating

CODE	DESCRIPTION	DIAMETER	LENGTH
94-1049SL-45	Slide Lock Handle	5 mm	45 cm



Snake Tooth (COBRA) Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
94-1069L	Direct Drive, Locking Ratchet on Demand	5 mm	33 cm



Stone Retrieval Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
95-1070	Direct Drive, Non-Locking	10 mm	33 cm

Symmetry Surgical Laparoscopic Instruments

360° rotation with flush port and Direct Drive system.
Laparoscopic instrumentation available in 45 cm length (-45).



Stone Spoon Retrieval & Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
94-1070	Direct Drive, Non-Locking	5 mm	33 cm



Scissors

Straight

CODE	DESCRIPTION	DIAMETER	LENGTH
94-2009	Direct Drive, Non-Locking	5 mm	33 cm



94-1051

Wave Atraumatic Grasping Forceps

Monopolar Handles

CODE	DESCRIPTION	DIAMETER	LENGTH
94-1051	Direct Drive, Non-Locking, Single-Action	5 mm	33 cm
94-1051-45	Direct Drive, Non-Locking, Single-Action	5 mm	45 cm
94-1051L	Direct Drive, Locking Ratchet on Demand, Single-Action	5 mm	33 cm


Symmetry Surgical Laparoscopic Instruments

360° rotation with flush port and Direct Drive system.
Laparoscopic instrumentation available in 45 cm length (-45).



Alligator Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
95-1052L	Direct Drive, Locking Ratchet on Demand	5 mm	33 cm
95-1052L-45	Direct Drive, Locking Ratchet on Demand	5 mm	45 cm



Allis Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
95-1078L	Direct Drive, Locking Ratchet on Demand	5 mm	33 cm



Babcock Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
95-1089L	Direct Drive, Locking Ratchet on Demand	5 mm	33 cm
95-1089L-45	Direct Drive, Locking Ratchet on Demand	5 mm	45 cm
95-1079L	Direct Drive, Locking Ratchet on Demand	10 mm	33 cm

Symmetry Surgical Laparoscopic Instruments

360° rotation with flush port and Direct Drive system.
Laparoscopic instrumentation available in 45 cm length (-45).



Claw Grasping Extractor Forceps

2x3 Heavy Teeth, Single-Action

CODE	DESCRIPTION	DIAMETER	LENGTH
95-1005L	Direct Drive, Locking Ratchet on Demand	5 mm	33 cm
95-1015L	Direct Drive, Locking Ratchet on Demand	10 mm	33 cm
95-1005L-45	Direct Drive, Locking Ratchet on Demand	5 mm	45 cm
95-1015L-45	Direct Drive, Locking Ratchet on Demand	10 mm	45 cm



DeBakey Grasping Forceps

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
95-1092L	Direct Drive, Locking Ratchet on Demand	5 mm	40 mm	33 cm
95-1092L-45	Direct Drive, Locking Ratchet on Demand	5 mm	40 mm	45 cm

Dorsey Grasping Forceps

Fenestrated

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
95-1031SL	Direct Drive, Slide Lock	5 mm	30 mm	34.5 cm

Symmetry Surgical Laparoscopic Instruments

360° rotation with flush port and Direct Drive system.
Laparoscopic instrumentation available in 45 cm length (-45).

Grasping Forceps

Blunt Tip




CODE	DESCRIPTION	DIAMETER	LENGTH
95-1032L	Direct Drive, Locking Ratchet on Demand	5 mm	33 cm



Hook Scissors

CODE	DESCRIPTION	DIAMETER	LENGTH
95-2020-45	Direct Drive, Non-Locking	5 mm	45 cm



Hunter Bowel Grasping Forceps

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
95-1093L	Direct Drive, Locking Ratchet on Demand	5 mm	20 mm	33 cm

Symmetry Surgical Laparoscopic Instruments

360° rotation with flush port and Direct Drive system.
Laparoscopic instrumentation available in 45 cm length (-45).



Maryland Dissecting Forceps

Curved

CODE	DESCRIPTION	DIAMETER	LENGTH
95-1067	Direct Drive, Non-Locking	5 mm	33 cm



Metzenbaum Scissors

Curved

CODE	DESCRIPTION	DIAMETER	LENGTH
95-2006	Direct Drive, Non-Locking	5 mm	33 cm
95-2006-45	Direct Drive, Non-Locking	5 mm	45 cm



Needle Nose Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
95-1045L	Direct Drive, Locking Ratchet on Demand	5 mm	32 cm

Symmetry Surgical Laparoscopic Instruments

360° rotation with flush port and Direct Drive system.
Laparoscopic instrumentation available in 45 cm length (-45).



Retraction Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
95-1056L	Direct Drive, Locking Ratchet on Demand, Single-Action	5 mm	33 cm



95-1048SL-DA-ROT

Retraction Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
95-1048	Single-Action	5 mm	33 cm
95-1048SL-ROT	Single-Action, Rotation	5 mm	35 cm
95-1048SL-ROT-45	Single-Action, Rotation	5 mm	45 cm
95-1048SL-DA	Double-Action	5 mm	35 cm
95-1048SL-DA-ROT	Double-Action, Rotation	5 mm	35 cm



95-1070

Stone Retrieval Forceps


CODE	DESCRIPTION	DIAMETER	LENGTH
95-1070	Direct Drive, Non-Locking	10 mm	33 cm

Symmetry Surgical Laparoscopic Instruments

360° rotation with flush port and Direct Drive system.
Laparoscopic instrumentation available in 45 cm length (-45).

Wave Grasping Forceps

Slide Lock Handle



CODE	DESCRIPTION	DIAMETER	LENGTH
95-1049SL		5 mm	33 cm
95-1049SL-45	Single-Action	5 mm	45 cm
95-1049SL-ROT	Single-Action, Rotation	5 mm	35 cm
95-1049SL-ROT-45	Single-Action, Rotation	5 mm	45 cm
95-1049SL-DA	Double-Action	5 mm	35 cm
95-1049SL-DA-45	Double-Action	5 mm	45 cm
95-1049SL-DA-ROT	Double-Action, Rotation	5 mm	35 cm

Wave Grasping Forceps



CODE	DESCRIPTION	DIAMETER	LENGTH
95-1051L	Direct Drive, Locking Ratchet on Demand, Single-Action	5 mm	33 cm

Wave Grasping Forceps

Double-Action



CODE	DESCRIPTION	DIAMETER	LENGTH
95-1053L	Direct Drive, Locking Ratchet on Demand	5 mm	33 cm
95-1053L-45	Direct Drive, Locking Ratchet on Demand	5 mm	45 cm

Symmetry Surgical Laparoscopic Instruments

360° rotation with flush port.


Laparoscopic instruments available in 45 cm length.



Tip

Atraumatic-Allis Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
97-1077RC	Ratchet Handle and Rotation	5 mm	35 cm



Tip

97-1079RC

Babcock Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
97-1089RC	Ratchet Handle and Rotation, Double-Action	5 mm	32 cm
97-1079RC	Ratchet Handle and Rotation, Double-Action	10 mm	37 cm



Tip

97-1030RC

Bottle Nose (Dolphin Nose) Dissecting Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
97-1030RC	Non-Locking Handle with Rotation, Fully-Insulated, Double-Action	5 mm	32 cm
97-1030RCL	Non-Locking Handle with Rotation, Fully-Insulated, Double-Action	5 mm	45 cm

Symmetry Surgical Laparoscopic Instruments

360° rotation with flush port.
Laparoscopic instruments available in 45 cm length.



Tip

Bowel Grasping Forceps

Fenestrated

CODE	DESCRIPTION	DIAMETER	LENGTH
97-1093RC	Ratchet Handle and Rotation	5 mm	32 cm



Cup Tooth Dissecting Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
97-1035RC	Non-Locking Handle with Rotation, Fully-Insulated, Double-Action	5 mm	32 cm



Tip

Cup Tooth Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
97-1055RC	Ratchet Handle and Rotation, Double-Action	5 mm	32 cm

Symmetry Surgical Laparoscopic Instruments

360° rotation with flush port.

Laparoscopic instruments available in 45 cm length.



Dorsey Grasping Forceps

Fenestrated

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
97-1031RC	Ratchet Handle and Rotation	5 mm	30 mm	32 cm

Duckbill Dissecting Forceps

Monopolar, Fully-Insulated



CODE	DESCRIPTION	DIAMETER	LENGTH
97-1040RC	Flattened Serrated Jaws, Non-Locking Handle with Rotation	5 mm	32 cm

Duval Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
97-1083RCL	Ratchet Handle and Rotation	10 mm	45 cm
97-1083RC	Ratchet Handle and Rotation	10 mm	32 cm

Symmetry Surgical Laparoscopic Instruments

360° rotation with flush port.
Laparoscopic instruments available in 45 cm length.




Hook Scissors

CODE	DESCRIPTION	DIAMETER	LENGTH
97-2020RC	Rotation, Fully-Insulated, Single-Action	5 mm	32 cm



Hunter Bowel Grasping Forceps

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
97-1090RC	Non-Locking Handle and Rotation, Fenestrated	5 mm	20 mm	32 cm
97-1090RCL	Non-Locking Handle and Rotation, Fenestrated	5 mm	20 mm	45 cm



Maryland Dissecting Forceps

Curved, Fully Insulated

CODE	DESCRIPTION	DIAMETER	LENGTH
97-1067RC	Non-Locking Handle and Rotation	5 mm	32 cm
97-1067RCL	Non-Locking Handle and Rotation	5 mm	45 cm

Symmetry Surgical Laparoscopic Instruments

360° rotation with flush port.

Laparoscopic instruments available in 45 cm length.

Maryland Tapered Dissecting Forceps

Monopolar

CODE	DESCRIPTION	DIAMETER	LENGTH
97-1063RC	Non-Locking Handle, Fully Insulated with Rotation	5 mm	35 cm



Metzenbaum Endo Scissors

CODE	DESCRIPTION	DIAMETER	LENGTH
97-2006RC	Non-Locking Handle with Rotation, Fully Insulated	5 mm	32 cm
97-2006RCL	Non-Locking Handle with Rotation, Fully Insulated	5 mm	45 cm



Mixer Dissecting Forceps

Double-Action

CODE	DESCRIPTION	DIAMETER	LENGTH
97-1066RC	Non-Locking Handle and Rotation, Fully Insulated	5 mm	33 cm
97-1068RC	Non-Locking Handle and Rotation, Fully Insulated	10 mm	33 cm
97-1068RCL	Non-Locking Handle and Rotation, Fully Insulated	10 mm	45 cm

Symmetry Surgical Laparoscopic Instruments

360° rotation with flush port.
Laparoscopic instruments available in 45 cm length.



Retraction Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
97-1056RC	Ratchet Handle and Rotation, Single-Action	5 mm	37 cm



Snake Tooth (COBRA) Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
97-1069RC	Ratchet Handle and Rotation	5 mm	37 cm



Swanstrom Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
97-1080RC	Ratchet Handle and Rotation	5 mm	37 cm


Symmetry Surgical Laparoscopic Instruments

360° rotation with flush port.
Laparoscopic instruments available in 45 cm length.




Tenaculum Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
97-1008RC	Ratchet Handle and Rotation	5 mm	33 cm



Wave Atraumatic Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
97-1051RC	Ratchet Handle and Rotation, Single-Action	5 mm	37 cm
97-1051DA-RC	Ratchet Handle and Rotation, Double-Action	5 mm	37 cm



Allis Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
97-1078RC	Ratchet Handle and Rotation	5 mm	32 cm

Symmetry Surgical RapidClean® Instruments

Three-piece modular system with choice of insulated, spring or ratchet handle; outer shaft and jaw assembly.



90-0030

Handle

For 90-Series Modular Instruments


CODE	DESCRIPTION
90-0030	Insulated with Cautery Post and Rotation Knob
90-0020	Non-Insulated with Spring Handle and Rotation Knob
90-0010	Non-Insulated with Ratchet and Rotation Knob



90-1115

Insulated Outer Shaft

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1115	With Flushport	5 mm	32 cm
90-1115L	With Flushport	5 mm	45 cm
90-1116	Without Flushport	5 mm	32 cm
90-1116L	Without Flushport	5 mm	45 cm
90-1117	Without Flushport	10 mm	32 cm
90-1117L	Without Flushport	10 mm	45 cm



90-1079B

Babcock Grasping Forceps

Double-Action

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1079R	Non-Insulated with Ratchet and Rotation Knob	10 mm	32 cm
90-1079B	Replacement Jaw Assembly		

Symmetry Surgical RapidClean® Instruments

Three-piece modular system with choice of insulated, spring or ratchet handle; outer shaft and jaw assembly.



Babcock Grasping Forceps
Double-Action

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1079RL	Non-Insulated with Ratchet and Rotation Knob	10 mm	45 cm
90-1079BL	Replacement Jaw Assembly		



Babcock Grasping Forceps
Double-Action

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1089R	Non-Insulated with Ratchet and Rotation Knob	5 mm	32 cm
90-1089B	Replacement Jaw Assembly		



Biopsy Forceps
Monopolar Single-Action

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1075R	Non-Insulated Ratchet Handle with Rotation Knob	5 mm	32 cm
90-1075B	Replacement Jaw Assembly		

Symmetry Surgical RapidClean® Instruments

Three-piece modular system with choice of insulated, spring or ratchet handle; outer shaft and jaw assembly.



Bottle Nose (Dolphin Nose)

Dissecting Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1030I	Monopolar, Double-Action, Insulated Handle	5 mm	32 cm
90-1030IL	Monopolar, Double-Action, Insulated Handle	5 mm	45 cm
90-1030B	Replacement Jaw Assembly		
90-1030BL	Replacement Jaw Assembly		



Claw Grasping Extractor Forceps

2x3 Heavy Teeth, Single-Action

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1015S	Spring Handle	10 mm	32 cm
90-1015B	Replacement Jaw Assembly		



Claw Grasping Extractor Forceps

2x3 Heavy Teeth, Single-Action

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1005S	Spring Handle	5 mm	32 cm
90-1005B	Replacement Jaw Assembly		

Symmetry Surgical RapidClean® Instruments

Three-piece modular system with choice of insulated, spring or ratchet handle; outer shaft and jaw assembly.



Cup Tooth Dissecting Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1035I	Insulated Handle, Double-Action	5 mm	32 cm
90-1035B	Replacement Jaw Assembly		



Duckbill Dissecting Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1040I	Monopolar, Fully-Insulated, Flattened Serrated Jaws, Double-Action with Insulated Handle	5 mm	32 cm
90-1040B	Replacement Jaw Assembly		



Fenestrated Dissecting Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1042IL	Insulated Handle with Rotation Knob	5 mm	45 cm
90-1042BL	Replacement Jaw Assembly		

Symmetry Surgical RapidClean® Instruments

Three-piece modular system with choice of insulated, spring or ratchet handle; outer shaft and jaw assembly.



90-1042B

Fenestrated Grasping Forceps

Double-Action

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1042R	Non-Insulated Ratchet Handle with Rotation Knob	5 mm	32 cm
90-1042B	Replacement Jaw Assembly		



90-1050B

Fundus Trimax Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1050R	Ratchet Handle, Double-Action	5 mm	33 cm
90-1050B	Replacement Jaw Assembly, Double-Action		



Tip

Maryland Dissecting Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1067I	Monopolar, Fully-Insulated, Double-Action with Insulated Handle	5 mm	32 cm
90-1067B	Replacement Jaw Assembly		

Symmetry Surgical RapidClean® Instruments

Three-piece modular system with choice of insulated, spring or ratchet handle; outer shaft and jaw assembly.



Tip

Maryland Dissecting Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1067IL	Monopolar, Fully-Insulated, Double-Action with Insulated Handle	5 mm	45 cm
90-1067BL	Replacement Jaw Assembly		



Tip

90-1065I

Maryland Dissecting Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1065I	Monopolar, Fully-Insulated, Double-Action with Insulated Handle	10 mm	32 cm
90-1065B	Replacement Jaw Assembly		



90-1066B

Tip

Mixer Dissecting Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1066I	Double-Action with Insulated Handle	5 mm	32 cm
90-1066B	Replacement Jaw Assembly		

Symmetry Surgical RapidClean® Instruments

Three-piece modular system with choice of insulated, spring or ratchet handle; outer shaft and jaw assembly.



Mixer Dissecting Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1068I	Double-Action with Insulated Handle	10 mm	32 cm
90-1068B	Replacement Jaw Assembly		



Snake Tooth (COBRA) Grasping Forceps

2 x 2 Teeth, Double-Action

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1069R	Non-Insulated Ratchet Handle with Rotation Knob	5 mm	32 cm
90-1069B	Replacement Jaw Assembly		



Spike Tooth Grasper

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1071R	Non-Insulated Ratchet Handle with Rotation Knob	5 mm	32 cm
90-1071B	Replacement Jaw Assembly		

Symmetry Surgical RapidClean® Instruments

Three-piece modular system with choice of insulated, spring or ratchet handle; outer shaft and jaw assembly.



Stone Retrieval Forceps

Single-Action, Ratchet Handle

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1070R	Spoon tips	10 mm	32 cm
90-1070B	Replacement Jaw Assembly		



Surgi-Shear 2-Piece Scissors

Double-Action, Curved Beak-Like Blades

CODE	DESCRIPTION	DIAMETER	LENGTH
90-2006I	Double-Action with Insulated Handle	5 mm	32 cm
90-2008B	Replacement Jaw Assembly with Outer Sheath		



Wave Atraumatic Grasping Forceps

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1051R	Ratchet and Rotation Knob	5 mm	32 cm
90-1051B	Replacement Jaw Assembly		

Symmetry Surgical RapidClean® Instruments

Three-piece modular design instruments with click and rotate.
Laparoscopic instruments available in 45 cm length.



Allis Atraumatic Grasping Forceps

Complete Instrument: Consists of Sheath, Replacement Jaw Assembly and Locking Carbon Fiber Handle

CODE	DIAMETER	LENGTH
96-1023R	5 mm	33 cm



Babcock Grasping Forceps

Complete Instrument: Consists of Sheath, Replacement Jaw Assembly and Locking Carbon Fiber Handle

CODE	DESCRIPTION	DIAMETER	LENGTH
96-1089R	DeBakey Teeth	5 mm	33 cm
96-1079R	DeBakey Teeth	10 mm	33 cm



Bottle Nose (Dolphin Nose) Grasping Forceps

Complete Instrument: Consists of Sheath, Replacement Jaw Assembly and Non-Locking Carbon Fiber Handle

CODE	DESCRIPTION	DIAMETER	LENGTH
96-1030NL	Monopolar, Double-Action	5 mm	33 cm

Symmetry Surgical RapidClean® Instruments


Three-piece modular design instruments with click and rotate.
Laparoscopic instruments available in 45 cm length.



Carbon Fiber Handle

For 96-Series Click and Rotating Modular Instruments

CODE	DESCRIPTION
96-0002NL	Non-Locking Carbon Fiber Handle
96-0003R	Locking Carbon Fiber Handle
96-0004R	Traditional Locking Ratchet Handle




96-1015Q

Claw Grasping Extractor Forceps

Complete Instrument: Consists of Sheath, Replacement Jaw Assembly and Locking Carbon Fiber Handle

CODE	DESCRIPTION	DIAMETER	LENGTH
96-1015R	2 x 3 Heavy Teeth, Single Action, Locking Carbon Fiber Handle	10 mm	33 cm



96-1035NL

Cup Tooth Dissecting Forceps

Complete Instrument: Consists of Sheath, Replacement Jaw Assembly and Non-Locking Carbon Fiber Handle

CODE	DESCRIPTION	DIAMETER	LENGTH
96-1035NL	Double-Action	5 mm	33 cm

Symmetry Surgical RapidClean® Instruments

Three-piece modular design instruments with click and rotate.
Laparoscopic instruments available in 45 cm length.



DeBakey Grasping Forceps

Complete Instrument: Consists of Sheath, Replacement Jaw Assembly and Locking Carbon Fiber Handle

CODE	DESCRIPTION	DIAMETER	TIP SIZE	LENGTH
96-1092NL	Double-Action	5 mm	40 mm	33 cm



Duckbill Dissecting Forceps

Complete Instrument: Consists of Sheath, Replacement Jaw Assembly and Non-Locking Carbon Fiber Handle

CODE	DESCRIPTION	DIAMETER	LENGTH
96-1022NL	Monopolar, Fully-Insulated Double-Action, with Insulated Handle, with Flattened Serrated Jaws	5 mm	33 cm




Fenestrated Dissecting Forceps

Complete Instrument: Consists of Sheath, Replacement Jaw Assembly and Non-Locking Carbon Fiber Handle

CODE	DESCRIPTION	DIAMETER	LENGTH
96-1042NL	Double-Action with Insulated Handle	5 mm	33 cm

Symmetry Surgical RapidClean® Instruments

Three-piece modular design instruments with click and rotate.
Laparoscopic instruments available in 45 cm length.



Hook Scissors

CODE	DESCRIPTION	DIAMETER	LENGTH
96-2020NL	Consists of Sheath, Replacement Jaw Assembly and Non-Locking Carbon Fiber Handle	5 mm	33 cm



Maryland Dissecting Forceps

Monopolar, Fully-Insulated, Double-Action, Non-Locking Carbon Fiber Handle

CODE	DESCRIPTION	DIAMETER	LENGTH
96-1067NL	Consists of Sheath, Replacement Jaw Assembly and Non-Locking Carbon Fiber Handle	5 mm	33 cm
90-1067BL	Replacement Jaw Assembly		



Mixer Dissecting

Complete Instrument: Consists of Sheath, Replacement Jaw Assembly and Non-Locking Carbon Fiber Handle

CODE	DIAMETER	LENGTH
96-1068NL	10 mm	33 cm

Symmetry Surgical RapidClean® Instruments

Three-piece modular design instruments with click and rotate.
Laparoscopic instruments available in 45 cm length.



Mixer Dissecting Forceps

Complete Instrument: Consists of Sheath, Replacement Jaw Assembly and Non-Locking Carbon Fiber Handle


CODE	DIAMETER	LENGTH
96-1066NL	5 mm	33 cm



Outer Sheath

For 96-Series Modular Instruments

CODE	DESCRIPTION	DIAMETER
96-0526	Rotation and Flush Port, for 33 cm Instruments	5 mm
96-0526-45	Rotation and Flush Port, for 45 cm Instruments	5 mm
96-0126	Rotation and Flush Port, for 10 cm Instruments	10 mm



Retraction Grasping Forceps

Complete Instrument: Consists of Sheath, Replacement Jaw Assembly and Locking Carbon Fiber Handle

CODE	DESCRIPTION	DIAMETER	LENGTH
96-1056R	Atraumatic Teeth	5 mm	33 cm

Symmetry Surgical RapidClean® Instruments


Three-piece modular design instruments with click and rotate.
Laparoscopic instruments available in 45 cm length.



Retraction Grasping Forceps

Complete Instrument: Consists of Sheath, Replacement Jaw Assembly and Locking Carbon Fiber Handle


CODE	DESCRIPTION	DIAMETER	LENGTH
96-1056R-45	Atraumatic Teeth	5 mm	45 cm



Snake Tooth (COBRA) Grasping Forceps

Complete Instrument: Consists of Sheath, Replacement Jaw Assembly and Locking Carbon Fiber Handle

CODE	DESCRIPTION	DIAMETER	LENGTH
96-1069R	2 x 2 Teeth, Double-Action	5 mm	33 cm




Stone Retrieval Forceps

Complete Instrument: Consists of Sheath, Replacement Jaw Assembly and Locking Carbon Fiber Handle

CODE	DESCRIPTION	DIAMETER	LENGTH
96-1070R	Single-Action	5 mm	33 cm

Symmetry Surgical RapidClean® Instruments

Three-piece modular design instruments with click and rotate.
Laparoscopic instruments available in 45 cm length.



Stone Spoon Retrieval & Grasping Forceps

Complete Instrument: Consists of Sheath, Replacement Jaw Assembly and Locking Carbon Fiber Handle

CODE	DESCRIPTION	DIAMETER	LENGTH
96-1072R	Ratchet Handle, Single-Action	10 mm	33 cm



Tenaculum Grasping Forceps

Complete Instrument: Consists of Sheath, Replacement Jaw Assembly and Locking Carbon Fiber Handle

CODE	DESCRIPTION	DIAMETER	LENGTH
96-1008R	Locking, 1 x 1 Sharp Teeth, Double-Action	5 mm	33 cm


Tenaculum Grasping Forceps

Complete Instrument: Consists of Sheath, Replacement Jaw Assembly and Locking Carbon Fiber Handle

CODE	DESCRIPTION	DIAMETER	LENGTH
96-1009R	Locking, 1 x 1 Sharp Teeth, Double-Action	10 mm	33 cm

Symmetry Surgical RapidClean® Instruments

Three-piece modular design instruments with click and rotate.
Laparoscopic instruments available in 45 cm length.




96-1053R

Wave Atraumatic Grasping Forceps

Complete Instrument: Consists of Sheath, Replacement Jaw Assembly and Locking Carbon Fiber Handle

CODE	DESCRIPTION	DIAMETER	LENGTH
96-1053R	Locking, Double-Action	5 mm	33 cm

Symmetry Ultra™ Kleppinger Type Bipolar Forceps



90-1018

Tip

Paddle Tong

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1018	Insert Only	5 mm	32 cm
90-1019	Insert Only	5 mm	45 cm

Symmetry Ultra™ Kleppinger Type Bipolar Forceps



Ring Handle

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1011	Ring Handle	5 mm	32 cm
90-1012	Ring Handle	5 mm	45 cm



Thumb Ring

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1013	Thumb Ring	5 mm	32 cm
90-1014	Thumb Ring	5 mm	45 cm



Waffle Tong

CODE	DESCRIPTION	DIAMETER	LENGTH
90-1016	Insert Only	5 mm	32 cm
90-1017	Insert Only	5 mm	45 cm

Electrodes

Laparoscopic Monopolar Electrodes

Single-Use, 10 per Box

CODE	DESCRIPTION	LENGTH	
90-9525-S	J-Hook	13 in	33 cm
90-9525-45S	J-Hook	17 ³ / ₄ in	45 cm
90-9522-S	L-Hook	13 in	33 cm
90-9522-45S	L-Hook	17 ³ / ₄ in	45 cm
90-9523-S	Spatula	13 in	33 cm
90-9523-45S	Spatula	17 ³ / ₄ in	45 cm



Laparoscopic Monopolar Electrode

Reusable

CODE	DESCRIPTION	LENGTH	
90-9525	J-Hook	12 ¹ / ₂ in	32 cm
90-9522	L-Hook	12 ¹ / ₂ in	32 cm
90-9523	Spatula	12 ¹ / ₂ in	32 cm
90-9526	Needle	12 ¹ / ₂ in	32 cm
90-9527	Button Tip	12 ¹ / ₂ in	32 cm



Laparoscopic Monopolar Suction/Irrigation Electrode

Reusable, 5 mm Shaft, with Valve, Fully Insulated

CODE	DESCRIPTION	LENGTH	
90-4005	J-Hook	12 ¹ / ₂ in	32 cm
90-4012	L-Hook	12 ¹ / ₂ in	32 cm
90-4015	Spatula	12 ¹ / ₂ in	32 cm
90-4010	Button Tip	12 ¹ / ₂ in	32 cm



Electrodes



Monopolar Cords

Single-Use, 10 per Box

CODE	DESCRIPTION	PLUG	LENGTH
40-2110	Endoscopic Instrument Connector, Straight	Valleylab-Style Generator	10 ft 3 m
40-2106	Endoscopic Instrument Connector, Angled 90°	Valleylab-Style Generator	10 ft 3 m



Monopolar Cords

Reusable, 1 Each

CODE	DESCRIPTION	PLUG	LENGTH
90018	Olsen® Threaded Instrument Connector	$\frac{5}{32}$ in (4 mm) Male Generator	10 ft 3 m
90020	Olsen® Threaded Instrument Connector	Valleylab-Style Male Generator	10 ft 3 m
90054	Laparoscopic Instrument Connector	$\frac{5}{32}$ in (4 mm) Male Generator	10 ft 3 m
90056	Laparoscopic Instrument Connector	Valleylab-Style Male Generator	10 ft 3 m
90065	3 mm Olympus-Style Snare Connector	$\frac{5}{32}$ in (4 mm) Male Generator	10 ft 3 m
90067	3 mm Olympus-Style Snare Connector	Valleylab-Style Male Generator	10 ft 3 m
50-2107	Endoscopic Instrument Connector, Straight	Valleylab-Style Male Generator	10 ft 3 m
50-2111	Endoscopic Instrument Connector, Angled 90°	Valleylab-Style Male Generator	10 ft 3 m

Sterilization Solutions



Prevac Base & Lid

CODE	DESCRIPTION	DIMENSIONS	
C22854	Endoscopic Length	22 x 8 ½ x 4 in	559 x 216 x 101 mm
C22855	Endoscopic Length	22 x 8 ½ x 5 in	559 x 216 x 127 mm
C22856	Endoscopic Length	22 x 8 ½ x 6 in	559 x 216 x 152 mm
C22858	Endoscopic Length	22 x 8 ½ x 8 in	559 x 216 x 203 mm
C24138	Full Length	24 x 12 x 8 in	610 x 305 x 203 mm



C13854-SP



C22856-S

Sterrad Base & Lid

CODE	DESCRIPTION	DIMENSIONS	
C22854-SP	Endoscopic Length	22 x 8 ½ x 4 in	558 x 215 x 101 mm
C22856-SP	Endoscopic Length	22 x 8 ½ x 6 in	558 x 215 x 152 mm

Sterilization Solutions




Ultra™ Container System

The lightweight Ultra™ Container System helps keep heavy instrument sets within AORN and AAMI limits. Designed for durability, it weighs 30% less than comparable sterilization systems.

KEY BENEFITS:

- Extruded rails reinforce the base and lid rims
- 13 inch (330 mm) width provides acceptable packaging for most loaner instrumentation trays
- Five sizes accommodate both general and specialty instruments




B2175-LAP

Large Laparoscopic Insert

CODE	DESCRIPTION	HOLDS	DIMENSIONS
B2175-LAP	Small Insert, For use with C22858	10 Instruments	21 x 7 ½ x 6 in 533 x 190 x 152 mm
B2312-LAP	Large Insert, For use with C24138	15 Instruments	23 x 12 x 6 in 584 x 305 x 152 mm



68-1599



68-1598

Laparoscopic Instrument Tray

CODE	DESCRIPTION
68-1599	Instrument Tray
68-1598	Pin Mat for Laparoscopic Instrument Tray