



Syringes

## SOL-M™ Syringe, 3-piece

- **High-transparency barrel** - makes it easy to see syringe content.
- **Silicone lubrication** - ensures smooth plunger movement.
- **Safe plunger backstop** - prevents accidental withdrawal of plunger.
- **Large and stable finger grip** - ensures stability and comfort during injection.
- **Luer Slip (concentric and eccentric), Luer Lock and Catheter options.**

### SOL-M™ Syringe without Needle - Blister



Capacity	Luer Lock	Luer Slip Concentric	Luer Slip Eccentric	Catheter
1 ml	180001** 180001PP*	180011*		
3 ml	180003*	180003ST*		
5 ml	180005*	180005ST*		
10 ml	180010*	180010ST*	180010ET*	
20 ml	180020*	180020ST*	180020ET*	
30 ml	180030*			
60 ml	180060*		180060ET*	180060CT***
100 ml				1800100CT***

\* Polypropylene \*\* Polycarbonate \*\*\* Polypropylene with Luer Slip adapter and tip cap

## SOL-M™ Syringe, 3-piece



Syringes



Capacity from 1 ml to 60 ml.

1 ml: polypropylene and polycarbonate options.

Syringes from 3 ml to 60 ml are made from polypropylene.

DEHP- and PVC- free, not made with natural rubber latex.

Sterile and single-use.

Sterilization method: ethylene oxide.