

Titanium Ligation System

Vclip[®] Ligating Clips & Clip Applicators

Reliable Ligation



PROVIDING SPECIALISTS WITH
INNOVATIVE
HEALTHCARE SOLUTIONS

WITH THE CENTRAL FOCUS ON
PATIENTS' SAFETY AND USERS' COMFORT




Vclip[®] ligating clips

Vclip[®] Ligating Clips are implantable medical devices. They are widely used in surgery for marking and ligating purposes. Medical practitioners value the product for its hemostatic features. The product offered by Grena is supplied with inner surreys that shall enhance its functions by assuring proper grip on the tissue and enabling its proper nourishment. Microstructural outer surface of clip legs and their divergence improves affixing of the clip in the applier. Vclip[®] Ligating Clips are manufactured of a medical grade titanium which ensures proper level of biocompatibility. The applied technology guarantees the product meets the requirements of even the most demanding clinical institutes.

Vclip[®] Ligating Clips are delivered in colour-coded cartridges to ensure easy identification. A cartridge includes 4, 6 or 10 clips that ensure economical and comfortable use. Product details necessary for its identification are also located on tracking labels. They can be placed in patient's history form. The product is disposable, EO sterilized.

KEY FEATURES:


- 6 sizes available: XS, S, SM, M, ML, L
- 4, 6 or 10 clips per cartridge
- Chevron-shaped
- Made of a medical grade titanium,
- Color-coded cartridges for easy identification
- Compatible with corresponding Vclip[®] Appliers
- Adhesive backing of the cartridge



Vclip[®]



REF	CLIP SIZE	COLOUR	COMPATIBLE WITH APPLIERS	PRODUCT	PCS PER CARTRIDGE/ BOX
0301-06XS	XS	○	50° - 60° V-shaped jaws cross-section	chevron-shaped titanium	6/30
0301-06S	S	●	50° - 60° V-shaped jaws cross-section	chevron-shaped titanium	6/30
0301-06S10	S	●	50° - 60° V-shaped jaws cross-section	chevron-shaped titanium	10/20
0301-06SM	SM	●	50° - 60° V-shaped jaws cross-section	chevron-shaped titanium	6/30
0301-06M	M	●	50° - 60° V-shaped jaws cross-section	chevron-shaped titanium	6/30
0301-06M10	M	●	50° - 60° V-shaped jaws cross-section	chevron-shaped titanium	10/20
0301-06ML	ML	●	50° - 60° V-shaped jaws cross-section	chevron-shaped titanium	6/20
0301-06ML04	ML	●	50° - 60° V-shaped jaws cross-section	chevron-shaped titanium	4/20
0301-06ML10	ML	●	50° - 60° V-shaped jaws cross-section	chevron-shaped titanium	10/15
0301-06L	L	●	50° - 60° V-shaped jaws cross-section	chevron-shaped titanium	6/20



10 CLIPS PER CARTRIDGE
Cost-effective alternative if you need more clips

4 CLIPS PER CARTRIDGE
Cost-effective alternative for simple procedures

Vclip[®] open surgery clip appliers

Vclip[®] open surgery appliers are available in a wide range of sizes and applier lengths. Made of the highest quality surgical stainless steel to meet the requirements of the most demanding professionals. Handle rings are color-coded matching cartridge hues.

KEY FEATURES:

- Available for XS, S, SM, M, ML and L non-absorbable Vclip[®] clips to be used on large spectrum of tissue sizes
- 3 lengths available: 15 mm, 20 mm and 28 mm
- Made of high quality stainless steel
- Colour-coded rings matching cartridge hue



REF	CLIP SIZE	COLOUR	TYPE	APPLIER LENGTH (CM)	PCS PER BOX
0301-07XS15*	XS	○	Open Surgery	15	1
0301-07XS20	XS	○	Open Surgery	20	1
0301-07XS28*	XS	○	Open Surgery	28	1
0301-07S15*	S	●	Open Surgery	15	1
0301-07S20	S	●	Open Surgery	20	1
0301-07S28*	S	●	Open Surgery	28	1
0301-07SM15*	SM	●	Open Surgery	15	1
0301-07SM20	SM	●	Open Surgery	20	1
0301-07SM28*	SM	●	Open Surgery	28	1
0301-07M15*	M	●	Open Surgery	15	1
0301-07M20	M	●	Open Surgery	20	1
0301-07M28*	M	●	Open Surgery	28	1
0301-07ML20	ML	●	Open Surgery	20	1
0301-07ML28	ML	●	Open Surgery	28	1
0301-07L20	L	●	Open Surgery	20	1
0301-07L28	L	●	Open Surgery	28	1

* Manufactured on request

Other lengths and jaw's angles available on request



Vclip[®] endoscopic clip appliers

Vclip[®] endoscopic clip appliers are available in a wide range of sizes. Made of the highest quality surgical stainless steel Grena appliers guarantees ergonomic and effective approach in endoscopic procedures. Handle rings are color-coded matching cartridge hues. Precise performance makes Grena appliers the leading choice for surgeons performing endoscopic procedures.

KEY FEATURES:

- Jaws activated by front handle for better stability,
- 360° shaft rotation by index finger allows one-hand operation
- Flushing port facilitates easy cleaning, thus eliminating potential contamination
- Made of high quality surgical stainless steel
- Non-detachable design for easier maintenance and durability



REF	CLIP SIZE	COLOUR	TYPE	WORKING LENGTH (CM)	TROCAR DIA-METER (MM)	PCS PER BOX
0301-07ME	M	●	Endoscopic	33	10	1
0301-07MLE	ML	●	Endoscopic	33	10	1
0301-07LE	L	●	Endoscopic	33	12	1

Omnifinger™ articulating endoscopic clip appliers

Omnifinger™ is an outstanding articulating endoscopic clip applier designed to deal with challenging cases. This innovative tool offers some unique features that can help surgeons to open the doors to new operating techniques.

KEY FEATURES:

- 60° of total angulation
- Articulation and 360° shaft rotation by index finger allows one-hand operation
- Flushing port facilitates easy cleaning, thus eliminating potential contamination
- Made of high quality surgical stainless steel



REF	CLIP SIZE	COLOUR	LOCK	WORKING LENGTH (CM)	TROCAR DIA-METER (MM)	PCS PER BOX
0301-07MEOMN	M	●	YES	36	10	1
0301-07MEOMNB	M	●	YES	47	10	1
0301-07MLEOMN	ML	●	YES	36	10	1
0301-07MLEOMNB	ML	●	YES	47	10	1

Vclip[®] endoscopic clip appliers 25°

Vclip[®] endoscopic applier, angled at 25 degrees, designed to provide precision and ease during surgical procedures. Our applier is made with high-quality materials to ensure durability and reliability in the operating room.

KEY FEATURES:

- Jaws activated by front handle for better stability
- 360° shaft rotation by index finger allows one-hand operation
- Flushing port facilitates easy cleaning, thus eliminating potential contamination
- Made of high quality surgical stainless steel
- Non-detachable design for easier maintenance and durability
- Colour-coded rings matching cartridge hue



REF	CLIP SIZE	COLOUR	TYPE	WORKING LENGTH (CM)	TROCAR DIA-METER (MM)	PCS PER BOX
0301-07MLEA25	ML	●	Endoscopic	33	10	1



Grena Ltd.

1000 Great West Road, Brentford,
London TW8 9HH, United Kingdom

EC REP MDML INTL LTD.

10 McCurtain Hill, Clonakilty, Co. Cork,
P85 K230, Republic of Ireland

CUSTOMER SERVICE CENTER NOTTINGHAM OFFICE, UK

☎ + 44 115 9704 800

☎ + 44 115 9704 800

✉ export@grena.co.uk

www.grena.co.uk

