

User Instructions

> Technical data

Name	Material
Weight	650g without handle / 800g with handle
Padding material	PU foam with PA hook and loop velour
Strap material	PA strap with PA hook and loop velour
Handle	latex-free foam material
Brace material	aluminium

> Size chart and article numbers

Name	Length of upper arm shell	Length of forearm shell	Circumference upper arm shell	Circumference forearm shell	Art.-No. left	Art.-No. right
Elbow Extension L/L	17 - 21 cm	20 - 24 cm	28 - 39 cm	22 - 30 cm	906LL-L	906LL-R
Elbow Extension L/L with handle	17 - 21 cm	20 - 24 cm	28 - 39 cm	22 - 30 cm	9061LL-L	9061LL-R
Elbow Extension M/M	17 - 20,5 cm	19,5 - 22 cm	20 - 30 cm	18 - 24 cm	906MM-L	906MM-R
Elbow Extension S/S	17 - 20 cm	17,5 - 20,5 cm	15 - 21 cm	14 - 19 cm	906SS-L	906SS-R
Elbow Extension Child Blue	17 - 20 cm	17,5 - 20,5 cm	15 - 21 cm	14 - 19 cm	906SS-BL	906SS-BR
Elbow Extension Child Orange	17 - 20 cm	17,5 - 20,5 cm	15 - 21 cm	14 - 19 cm	906SS-OL	906SS-OR
Elbow Extension Child Pink	17 - 20 cm	17,5 - 20,5 cm	15 - 21 cm	14 - 19 cm	906SS-PL	906SS-PR
Elbow Extension Child Green	17 - 20 cm	17,5 - 20,5 cm	15 - 21 cm	14 - 19 cm	906SS-GL	906SS-GR

Name	Art.-No. left	Art.-No. right
Rotatable handle	993A-HL	993A-HR