

User Instructions

> Technical data

Name	Material
Weight	1048g
Padding material	PU foam with PA hook and loop velour
Strap material	PA strap with PA hook and loop velour
Brace material	aluminium

> Size chart and article numbers

Name	Length of thigh shell	Length of lower leg shell	Circumference of thigh	Circumference of lower leg	Art.-No. left	Art.-No. right
Knee Extension L/L	23 - 25,5 cm	26,5 - 29 cm	36 - 66 cm	32 - 47 cm	901LL-L	901LL-R
Knee Extension L/M	23 - 25,5 cm	26,5 - 29 cm	33-43 cm	24 - 35 cm	901LM-L	901LM-R
Knee Extension M/L	18,5 - 21 cm	20,5 - 23,5 cm	36 - 66 cm	32 - 47 cm	901ML-L	901ML-R
Knee Extension S/S	17,5 - 20 cm	21,5 - 25 cm	30 - 40 cm	22 - 32 cm	901SS-L	901SS-R
Knee Extension XS/XS	16,5 - 18,5 cm	18,5 - 21,5 cm	22 - 32 cm	20 - 25 cm	901XSXS-L	901XSXS-R
Knee Extension Child Blue	17,5 - 20 cm	21,5 - 25 cm	30 - 40 cm	22 - 32 cm	901SS-BL	901SS-BR
Knee Extension Child Orange	17,5 - 20 cm	21,5 - 25 cm	30 - 40 cm	22 - 32 cm	901SS-OL	901SS-OR
Knee Extension Child Pink	17,5 - 20 cm	21,5 - 25 cm	30 - 40 cm	22 - 32 cm	901SS-PL	901SS-PR
Knee Extension Child Green	17,5 - 20 cm	21,5 - 25 cm	30 - 40 cm	22 - 32 cm	901SS-GL	901SS-GR
Knee Extension Infant Blue	16,5 - 18,5 cm	18,5 - 21,5 cm	22 - 32 cm	20 - 25 cm	901XSXS-BL	901XSXS-BR
Knee Extension Infant Orange	16,5 - 18,5 cm	18,5 - 21,5 cm	22 - 32 cm	20 - 25 cm	901XSXS-OL	901XSXS-OR
Knee Extension Infant Pink	16,5 - 18,5 cm	18,5 - 21,5 cm	22 - 32 cm	20 - 25 cm	901XSXS-PL	901XSXS-PR
Knee Extension Infant Green	16,5 - 18,5 cm	18,5 - 21,5 cm	22 - 32 cm	20 - 25 cm	901XSXS-GL	901XSXS-GR