

User Instructions

> Technical data

Name	Material
Weight	1264g / Child Brace 730g
Padding material	PU foam with PA hook and loop velour
Strap material	PA strap with PA hook and loop velour
Brace material	aluminium
Sole	EVA foam

> Size chart and article numbers

Name	Length of thigh shell	Circumference of the calf	EU- shoe sizes	Sole Length	Art.-No. left	Art.-No. right
Ankle L/L	36 - 38,5 cm	32 - 47 cm	44 - 47	30 cm	920LL-L	920LL-R
Ankle L/M	36 - 38,5 cm	32 - 47 cm	40 - 43	27,5 cm	920LM-L	920LM-R
Ankle L/S	36 - 38,5 cm	32 - 47 cm	37 - 39	26 cm	920LS-L	920LS-R
Ankle S/SX	27 - 31 cm	22 - 32 cm	24 - 34	23 cm	920SXS-L	920SXS-R
Ankle Child Blue	27 - 31 cm	22 - 32 cm	24 - 34	23 cm*	920SXS-BL	920SXS-BR
Ankle Child Orange	27 - 31 cm	22 - 32 cm	24 - 34	23 cm*	920SXS-OL	920SXS-OR
Ankle Child Pink	27 - 31 cm	22 - 32 cm	24 - 34	23 cm*	920SXS-PL	920SXS-PR
Ankle Child Green	27 - 31 cm	22 - 32 cm	24 - 34	23 cm*	920SXS-GL	920SXS-GR

- * The foot bedding (insole) has to be made individually by the specialised trade in order to guarantee an optimal bedding of the child's foot. The special construction of the sole part of the paediatric CDS® Ankle Brace facilitates an individual adjustment through the orthopaedic technician.